


Agilent OpenLab CDS Quick Reference Sheet



www.agilent.com/chem/openlabcds

Edition 04/2022
Document Number D0002171en Rev. B

OpenLab CDS – Basic Information

Terms	
ACQUISITION METHOD	Includes the instrument settings for an experiment. Edited in Acquisition .
AUDIT TRAIL	A representation of changes to a record that includes 1) identity of the user who made the change, 2) date/time of the change, 3) description of the change, and 4) reason for the change. Configure audit trails in OpenLab Control Panel .
LAYOUTS	Defines how and which information is displayed. Use preset Layouts or customize your personal layout. Changes in layouts are saved per user. Use Reset function to return to default layouts.
LINKED METHOD	Processing method that is assigned to an injection.
PROCESSING METHOD	Contains the information and parameters needed to process the data and generate results. Edited in Data Analysis .
PROJECT	Use to organize your data, e.g. by instrument, laboratory, or study. Apply individual access rights to each project
REPORT TEMPLATE	Defines the layout of a report. Edited in Data Analysis .
RESULT SET	Collection of raw data, methods, and injection list. Create custom result sets by combining single samples or sequence injections, and reprocess together.
SEQUENCE CREATION TEMPLATE	Defines a flexible pattern for creating new sequences. Useful for creating cyclic sequence. Edited in Acquisition .
SEQUENCE TABLE	Identifies the order of runs to be included in a sequence. Edited in Acquisition .

Help and Learning

GETTING STARTED - OPENLAB HELP AND LEARNING

The interactive **e-Introduction** modules provide step-by-step instructions on basic workflows and help you to get familiar with OpenLab CDS. Access the modules from the desktop by clicking one of the **OpenLab Help and Learning** icons.


Use **F1** to get **context sensitive help**. Check Help and Learning Online for the latest localized versions of manuals and modules at <https://openlab.help.agilent.com>



AGILENT COMMUNITY

<https://community.agilent.com>
Collaborate with others about applications, discuss Agilent products, and find in-depth documents and videos relevant to Agilent.

MENU COMMANDS



Find commands in the ribbon. **Data Analysis**: contextual ribbon tabs (highlighted in blue) offer commands for the selected window.

Prepare, Install, Configure

Prepare:
Run the System Preparation tool, and check the Requirements guide for hardware, software and network requirements.

Install, configure:
Click **setup.exe** on the USB stick to run the installation wizard.

OpenLab CDS – Quick Reference Sheet

Access Instruments and Projects



Acquire Data



Analyze and Report Data



CONTROL PANEL



Use the Desktop icon to launch OpenLab Control Panel.

MANAGE USERS

ROLE: SYSTEM ADMIN

- Assign users to groups
- Assign roles to groups
- Use default roles or assign privileges to roles (use role **Chemist** for Acquisition and Data Analysis)
- Edit security policy
- Manage licenses

PROJECTS

ROLE: PROJECT ADMIN

Create projects to organize your data:

- Activate method audit trails
 - Apply project-specific options
 - Add custom parameters for samples and compounds
- Create shortcut for Data Analysis with a project.
- Launch Data Analysis for a project.



INSTRUMENTS

ROLE: INSTRUMENT ADMIN

Create and configure instruments
Monitor instruments



Create shortcut for acquisition with an instrument.

Launch Acquisition for an instrument.



CREATE ACQUISITION METHOD

Define acquisition and instrument parameters for runs using this method.

Download to instrument to use new parameters immediately.



RUN A SINGLE SAMPLE

- Enter sample parameters
- Include acquisition method and (optional) processing method
- Run



RUN A SEQUENCE

- Create a new sequence
- Open and edit a previously saved sequence
- Import sequence from a result set or a CSV file (also possible via drag and drop)
- Click **Apply Template** to apply a sequence creation template



SEQUENCE

Edit sequence table:

Right-click and use **Fill-down** to edit rows, **copy/paste** entire rows, **pin columns** at the front

Edit injections and sequence parameters:

- For automated processing and reporting: enter existing processing method
- Define name of the result set

SEQUENCE CREATION TEMPLATE

Create a new template in the **Sequence Creation Template** window:

Sequence Creation Template

Use for cyclic sequences with varying numbers of samples.



STATUS

View Run Queue, Instrument status, Online signals



LAUNCH DATA ANALYSIS

Click icon in the **History** tab of the **Run Queue** to review completed single injections or sequences in Data Analysis.

REVIEW DATA

- To load data, double-click node in the **Data Selection** view
- Pin injections or result sets for comparison
- Use **Peak Explorer** to view a sequence at a glance and identify trends or artifacts.
- Create custom result sets based on existing injections (right-click nodes in the injection tree to add them, click **Create new Result Set** in the ribbon).



Use the **Integration Optimizer** to quickly adjust peak baselines.

INJECTIONS AND PROCESSING METHODS



A processing method may already be linked (assigned) to the data.

If not, create new master method

Link injections to the master processing method



A result set method (stored in result set folder) is generated

on linking the master method to the injections.

Reprocess (automatic after linking, or manually)



EDIT PROCESSING METHOD

Processing methods are divided in sections: e.g.

- **Integration Events** for all or for specific signals
- **Compounds:**
Add compounds: Right-click a peak, or right-click the compound table.
- **Calibration:** Set number of calibration levels in the **General** tab
- **Spectra:** MS or UV reference spectra for confirmation or purity checks
- **Extraction > Spectrum:** Background correction settings
- **Tools > Custom Calculation:** Extend default calculations

REPORT AND PRINT



Import default templates into your project

(start customizing templates from here):

Data Selection > Import Default Templates



Print results manually:

In Data Analysis, edit the processing method **Reports** section, and link the method. Click **Print all** in the ribbon.

Print processing method:

Go to **Data Processing > File** ribbon tab > **Print Method**

Print individual PDF:

Go to **Reporting > Preview > Save as PDF**



Sequence summary report: submenu **Save result set summary report**