

6


Troubleshooting


6.1 Fault Diagnosis and Corrective Actions

This chapter describes the causes of problems and the corrective action to take when they occur. For more details on the action to take, refer to the reference pages cited in the tables below.

If the problem is not resolved after taking the corrective action described here, or if the fault is considered to be one not described here, contact your Shimadzu representative.

■ General

Problem	Cause / Corrective Action
<p>The LEDs do not light when the power is turned ON.</p>	<p>Check if the power cable is correctly connected and if single-phase 230 V power is being supplied to the distribution panel. (See "8.1.2 Connecting the Instrument" P.189.)</p> <p>Status immediately after turning the power ON</p>  <p>The diagram shows a line drawing of the instrument's front panel. A callout box highlights five LEDs labeled HV, GAS, HEATER, STATUS, and POWER. The STATUS LED is shown with radiating lines, indicating it is lit. Labels 'Lit in green' and 'Flashing in green' point to the STATUS LED.</p>
<p>The instrument does not connect.</p>	<p>Check the connection of the USB cable. (See "2.2 Names of Parts" P.14.)</p> <p>After checking that the power to the instrument (MS, LC) is ON, restart the PC.</p> <p>Check the environment settings in LabSolutions.</p>

Problem	Cause / Corrective Action
<p>The detector and heater cannot be turned ON. If they are set to ON they immediately turn back OFF again.</p>	<p>The vacuum system is not operating normally. If the vacuum system is not in the "ready" status the MS detector and heater cannot be turned ON.</p> <p>An interlock has been actuated. Interlock functions operate as a measure to prevent danger when people are working. If either of the conditions indicated in the figure below applies, the MS detector (high voltage) and heater are turned OFF.</p> 
<p>The instrument stops and the "STATUS" LED lights in red.</p>	<p>When an abnormality of the vacuum system has been detected, the instrument is automatically stopped by a protective function. The probable causes are as follows.</p> <ul style="list-style-type: none"> • A large vacuum leak • An abnormality of the vacuum pump (rotary pump, triple inlet turbo molecular pump) • An abnormality of the leak valve <p>If this happens, temporarily exit LabSolutions and turn OFF the power to the PC and MS. Checking for vacuum leaks and, after effecting repairs (see "7.12 Vacuum Leak Check" P.176), turn the power to the MS back ON, then turn the power to the PC ON and start up LabSolutions.</p>

Problem	Cause / Corrective Action
The "STATUS" LED is flashing in green.	<p>When the triple inlet turbo molecular pump is still in preparation, or the degree of vacuum has not reached the target value, the LED flashes in green.</p> <p>The LED flashes in green when in standby mode.</p> <p> NOTE When starting and stopping the instrument, except when it is started other than by an automatic start (i.e., manually), the ion gauge has to be set to ON.</p> <p> NOTE When the instrument has been stopped for a long time it takes some time to evacuate it and the "Ready" status may not be reached within the prescribed time. In this case the "STATUS" LED on the front of the instrument will continue to flash in green, but this does not indicate an abnormality. When the pressure has dropped the "Ready" status will be attained and the "STATUS" LED will light in green.</p> <p>If the ion gauge cannot be turned ON, or if the ion gauge is turned on but the "STATUS" LED continues to flash in green for longer than three hours, there could be a vacuum leak.</p> <p>Take appropriate action by referring to "7.12 Vacuum Leak Check" P.176.</p>
The "STATUS" LED is flashing in green and the LabSolutions screen indicates standby status.	The DL is not mounted. Mount the DL.
	The DL connector is disconnected. See "7.9 Replacing the DL" P.163 and check that the DL connector is connected properly.
	There is a broken wire. See "7.9 Replacing the DL" P.163 and replace the DL.
	Check the Pirani gauge value on the instrument status. If the value is more than 300 Pa, perform the procedure described in "7.12 Vacuum Leak Check" P.176 .
The buzzer is sounding.	Check the cause of the error, which will be displayed on the LabSolutions screen.
The buzzer sounds in the pattern "pip pipee pip pipee" (sounding twice per second).	<p>A liquid leakage has been detected at the leak tray inside the instrument's front door.</p> <p>When a liquid leakage error occurs, analysis is stopped. Eliminate the cause of the liquid leakage.</p>
The buzzer sounds in the pattern "pip pip pip pip" (sounding four times per second).	<p>This indicates overheating of the DL, APCI heater and heated block. The heater is turned OFF automatically. This indicates that, as the result of some abnormality, the temperature of any part of the DL, heated block or APCI heater, has risen to the upper limit value for heat resistance. When the upper limit temperature is reached the heaters are automatically turned OFF; once the temperature has fallen below the upper limit value the buzzer stops sounding. Once the temperature has dropped below 50 °C, check that the DL and APCI connectors are attached properly.</p>
	Overheating and liquid leakage are detected at the same time.
The message "Ion gauge error" is displayed.	<p>The filament of the ion gauge vacuum meter has broken.</p> <p>Contact your Shimadzu representative if the message appears even after turning the ion gauge ON again.</p>
The message "Pirani gauge error" is displayed.	<p>The filament of the Pirani gauge vacuum meter has broken.</p> <p>Contact your Shimadzu representative if the message appears even after checking the connection of the Pirani gauge again.</p>
The source window cannot be closed.	<p>The O-ring of the source window is lifted.</p> <p>Correctly fit the O-ring into the groove.</p>

■ Data Acquisition

Problem	Cause / Corrective Action	See
The ion intensity is unstable or low.	ESI probe	
	The capillary ASSY is blocked.	"7.2.1 Replacing the Tubing" P.122
	The tip of the capillary ASSY is projecting too far beyond the ESI probe, or is too recessed.	"7.2.1 Replacing the Tubing" P.122
	The ESI probe is substantially out of place.	"3.5.1 Mounting the ESI Probe" P.55
	The probe current (interface current) is too high. (It can become unstable as the result of discharge.)	Lower the interface voltage.
	The high-voltage cable is not connected.	"3.5.1 Mounting the ESI Probe" P.55
	The high voltage is not being supplied properly (broken wire).	Check the analysis method parameters and the tuning file.
	APCI probe	
	The APCI pipe ASSY is blocked.	"7.3.1 Replacing the APCI Pipe" P.129
	The corona needle is substantially displaced from the center.	"3.6.1 Mounting the Corona Needle for APCI" P.60
	The corona needle is soiled.	"7.3.3 Cleaning and Replacement of the APCI Corona Needle" P.136
	The probe current (interface current) is too high. (It can become unstable as the result of discharge.)	Lower the interface voltage.
	The probe has not been heated properly.	"3.6.2 Mounting the APCI Probe" P.63
	The heater cable and/or high-voltage cable are not connected.	
Broken wire in the high-voltage cable, heater cable or sensor.		

Problem	Cause / Corrective Action	See
The ion intensity is unstable or low.	DL	
	The DL is blocked (the PG value is 50 Pa or lower).	"7.9 Replacing the DL" P.163
	The DL is soiled.	"7.9 Replacing the DL" P.163
	The probe has not been heated properly.	Check the analysis method parameters and the tuning file. If the situation does not improve, see "7.9 Replacing the DL" P.163.
	The wiring has not been done properly. (Parts are shorted.)	"7.9 Replacing the DL" P.163
	Detector	
	The detector voltage is too low.	
	The detector has deteriorated.	
	The attenuation is set low.	Adjust the attenuation. "Attenuation settings" P.113
	Heated block	
	The heated block is soiled.	"7.7 Cleaning the Spray Unit" P.159
	Lens system	
	The skimmer is soiled.	"7.5 Maintenance of the Lens System" P.142
	The multipole is soiled.	"7.5 Maintenance of the Lens System" P.142
	The entrance lens is soiled.	"7.5 Maintenance of the Lens System" P.142
The skimmer or the entrance lens is blocked.	"7.5 Maintenance of the Lens System" P.142	
The base level (background level) is too high.	Soiling of the DL	"7.9 Replacing the DL" P.163
	Soiling of the heated block	"7.7 Cleaning the Spray Unit" P.159
	Soiling of the spray unit	"7.7 Cleaning the Spray Unit" P.159
	Soiling of the LC	Clean the mobile phase.
	Soiling of the ESI and/or APCI probe	Clean the mobile phase.
	Soiling due to the standard sample for auto tuning	"7.7 Cleaning the Spray Unit" P.159
	Soiling of the mobile phase	Replace the mobile phase.
	Influence of background noise originating from the mobile phase	
	Soiling of the tubing	Clean the flow line, or replace the tubing.

Problem	Cause / Corrective Action	See
The peaks are too wide.	There is dead volume in the tubing section (section that connects to the ionization probe).	Repair the connection.
	There is dead volume in the connecting section downstream of the injector.	Repair the connection.
	The tubing bore is too large (make the internal diameter downstream of the injector 0.13 mm).	Replace the tubing.
	The cut face of the tubing is inclined.	Replace the tubing.
	Soiling of the injector	Rinse the injector.
	The ESI probe is out of position.	"3.5.1 Mounting the ESI Probe" P.55
The baseline undulates.	Air is contained in the tubing.	Purge the LC pump.
	Air is contained in the mobile phase.	Degas the mobile phase.
	The temperature of the DL, heated block or APCI heater is fluctuating.	
The pump pressure is too high.	Blockage of the ESI capillary. Blockage of the APCI or SUS pipe	"7.2.1 Replacing the Tubing" P.122 "7.3.1 Replacing the APCI Pipe" P.129
	Blockage of the tubing	"7.2.1 Replacing the Tubing" P.122 "7.3.1 Replacing the APCI Pipe" P.129
	Blockage at the LC side	Check the LC side and clear the blockage.

Problem	Cause / Corrective Action	See
Auto tuning has failed.	The wrong standard sample was used.	"8.3.1 Method for Preparing the Standard Sample" P.213
	The standard sample is old (the spectrum has changed).	"8.3.1 Method for Preparing the Standard Sample" P.213
	There is a leak in the tubing of the standard sample introduction unit.	"8.1.5 Installing the Standard Sample" P.198
	The resistance tube is clogged.	"7.6.2 Replacing the Resistance Tube" P.156
	There is not a sufficient volume of standard sample in the standard sample bottle.	"7.6.1 Replacing the Standard Sample" P.154
	The ESI probe is out of position.	Adjust the ESI probe position to gain a stronger target m/z intensity. In particular, check that negative ions are m/z 1007.
	An interference peak exists around the target m/z.	After washing the ESI probe tubing with water or methanol, check that there is no interference peak around the target m/z in the manual tuning window. "8.2.1 Checking Ions in the [MS Tuning] Window" P.201
	The target m/z intensity decreases when using substances such as TFA.	"Interpreting tuning results" P.90 in "3.8.2 Starting Auto Tuning"
The peaks are split.	The ion intensity has exceeded 20 million. (The signal strength is too high, causing saturation.)	The sample concentration is too high. Dilute the sample before use.
	There is dead volume in the tubing.	Repair the tubing connection.
	There is a problem with the column.	Replace the column.
The peak top is flat.	The ion intensity has exceeded 20 million. (The signal strength is too high, causing saturation.)	Adjust the attenuation. "Attenuation settings" P.113

Problem / Cause	Corrective Action / See
Software	
The tuning file is not appropriate. An old tuning file is being used.	"3.8.2 Starting Auto Tuning" P.82
The method file is not appropriate.	Refer to the LabSolutions Getting Started Guide or Operators Guide.
The batch file is not appropriate.	
The lens setting conditions are wrong.	"3.8.2 Starting Auto Tuning" P.82
Tubing	
There is a liquid leakage. (The pump pressure is too low.)	Mount the ionization probe. Identify and retighten any locations at which liquid is leaking.
There is a gas leak.	"8.1.3 About the Gas Used" P.193
	"3.5.1 Mounting the ESI Probe" P.55 "3.6.2 Mounting the APCI Probe" P.63
The pumping is unstable.	Purge the LC pump and check the check valve (see the LC pump instruction manual).
The gas is not flowing.	"8.1.3 About the Gas Used" P.193
The mobile phase is not flowing.	

■ Attenuation settings


If signal saturation occurs due to excessive ion intensity, attenuation adjustment may solve the problem. If the problem persists, dilute the sample.

- 1 Click [Attenuation].
The [Attenuation] window is displayed.

The screenshot shows the Realtime Analysis software interface. The main window displays three chromatograms (LC, MS, and ALL) and the Instrument Parameters View. The Attenuation window is highlighted in the Instrument Parameters View, showing the Conversion Dynode voltage set to 6 kV.

Item	Value	Setting	Units
Nebulizing Gas Flow	---	3.0	L/min
Drying Gas Flow	---	15.0	L/min
Interface	DUIS - ESI		
Interface Voltage		0.0	kV
Interface Current	0.1		uA
DL Temperature	186	250	C
Heat Sink Temperature	227	400	C
Detector Voltage		0.00	kV
IG Vacuum	1.8e-003		Pa
PG Vacuum	1.3e+002		Pa
CID Gas	230	230	kPa
Mode	Isocratic Flow	Isocratic Flow	
Total Flow	0.0000	0.0000	mL/min
Flow Control	2.0	2.0	Hz
Pump A Flow	0.0000	0.0000	mL/min
Pump B Flow	0.0000	0.0000	mL/min
Pump A Pressure	0.0		MPa
Pump B Pressure	0.0		MPa
Pump A Degasser	Not Connected		kPa
Pump B Degasser	Not Connected		kPa
Oven Temperature	31.0	40	C
Temperature Limit(Maximum)	65	65	C
Injection Volume			uL

- 2 Adjust the conversion dynode voltage.
Set a voltage that will reduce the signal intensity to below 20 million. The default value is 6 kV.



6.2 Operation Flow for Automatically Starting the Vacuum System




Fig. 6-1

**NOTE**

When the instrument has been stopped for a long time it takes some time to evacuate it and the "Ready" status may not be reached within the prescribed time. In this case the "STATUS" LED on the front of the instrument will continue to flash in green, but this does not indicate an abnormality. When the pressure has dropped the "Ready" status will be attained and the "STATUS" LED will light in green.

6.3 Operation Flow for Automatically Stopping the Vacuum System




Fig. 6-2

**NOTE**

To turn the vacuum system off manually, turn each of the heaters OFF and, when the temperature has fallen to 100 °C or lower, stop the vacuum pump.

6.4 Power Outages


If the 230 V AC power supply to the LCMS-8030/LCMS-8040 has been cut off and LabSolutions is set to the automatic restart mode (i.e. if [Vacuum Restart Mode] is checked among the recovery modes in the [System Control] window), the vacuum system will automatically start when the power is restored.



NOTE

To set the automatic recovery mode, select [System Control] from the assistant bar in the LabSolutions analysis window and check the [Vacuum Restart Mode] box in the [System Control] window.

If automatic restart is not set, start the instrument by referring to "[3.1 Starting the Instrument](#)" P.41.



7

Maintenance

7.1 Periodic Inspections and Servicing

In order to use this product safely, carry out periodic inspections.

Conduct everyday maintenance and management of performance as appropriate.

For details on maintenance inspection contracts, consult your Shimadzu representative.



WARNING



Instructions

Before starting maintenance or servicing work, turn the power to the product OFF.

Otherwise there will be a risk of fire, electric shock and instrument failure.



CAUTION



Prohibitions

Never remove the main body cover.

This may cause injury or damage the instrument.

If repair work that necessitates removal of the main body cover is required, consult your Shimadzu representative.



Instructions

- When replacing parts, use the parts specified by Shimadzu.
Using parts other than those specified could lead to breakage of parts, injuries or instrument failure.
- When carrying out periodic inspections, observe the cautions on each operation.
- Wear clean gloves when carrying out maintenance work on the ion source unit, the lens system or the interface.
- Be sure to wipe off any soiling on the tools used for maintenance work on the ion source unit, lens system or interface using gauze impregnated with acetone.
Using soiled tools or gloves will cause background noise.
- Place removed parts on clean fabric such as gauze.
If you place parts on soiled fabric they will become soiled.

7.1.1 List of Periodic Inspections and Servicing

 CAUTION	
	The replacement and inspection frequencies indicated in this list are merely for guidance. They are not guaranteed intervals.
Instructions	There is some latitude in the appropriate frequency depending on the conditions of use.

■ Replacement

Inspection/Service Operation or Part	Part No.	Replacement/Inspection Frequency				Remarks
		Every 6 Months	Every Year	Every 2 Years	Every 3 Years	
ESI probe  P.122						
Capillary tube body	S225-14915		○			
Capillary ASSY	S225-14948-91		○			
PEEK tube	S228-32999-01		○			
Elbow, KQ2L23-M5	S035-60724-21				○*	
O-ring, 4D-S8	S036-19004-05				○*	Capillary holder
APCI probe  P.128						
Heater unit ASSY	S225-15619-41		○			
Corona needle (for both APCI and DUIS)	S225-15877-92				○	
APCI pipe ASSY	S225-15845-91		○			
Half union, KQ2H01-M5	S035-60725-01				○	
DUIS  P.140						
Corona needle (for both APCI and DUIS)	S225-15877-92				○	
Interface  P.162, P.163, P.166						
DL ASSY	S225-15718-91	○				
Orifice	S225-15479				○*	
Sampling cone	S225-15487				○*	
M6PEEK washer for heater flange	S023-65106-01			○		
O-ring, 4D-S100	S036-19004-53				○*	Heater flange
O-ring, 4D G145	S036-12527				○*	Front door
Lens system  P.177						
Multipole 1	-					
Multipole 2	-					
O-ring, AS568A-253	S036-15552-53				○*	Door unit

Inspection/Service Operation or Part	Part No.	Replacement/Inspection Frequency				Remarks
		Every 6 Months	Every Year	Every 2 Years	Every 3 Years	
Vacuum system						 P.246, P.247
IG GAUGE	S225-09490-01				○*	
Gasket 14	S261-00207-02				○*	IG unit
Pirani gauge filament	S225-20310-91				○*	
RP oil #46 (4L)	S017-30163-02		○			Changed 3 times per year 1.5 L used in each oil change
SI unit ASSY						 P.154
PEEK frit	S228-48607-91				○*	
Resistance tube, alone	S225-15873-91			○		
Spacer, FKM	S225-15697-01		○			
Detector						 P.31
EM, MS644	S225-14168-01				○**	
N ₂ gas related						 P.178
FEP tube for nebulizer gas, 1/16	S225-14255-41			○		Cut off 1 cm from the end of the tube once every 2 years.
Sample						 P.213
Sample for auto tuning (200 mL)	S225-14122-01		○			

**NOTE**

The asterisks in the "Replacement/Inspection Frequency" column indicate that the replacement parts are included in the PM 3-year maintenance part kit (P/N S225-15869-42). Work indicated with ** must be performed by field engineers.

Overhaul

Inspection/Service Operation or Part	Part No.	Replacement/Inspection Frequency				Remarks
		Every 6 Months	Every Year	Every 2 Years	Every 3 Years	
Vacuum system						 P.245
RP E2M28 230V	-			○		***
Triple inlet turbo molecular pump body and SPLIT FLOW310 set	-				○	

**CAUTION****Instructions**

Perform overhauls.

If a rotary pump for which "***" is indicated under "Overhaul" is not overhauled, it may cause an accident or performance may decrease, resulting in damage to the instrument.

■ Cleaning

Inspection/Service Operation or Part	Part No.	Replacement/Inspection Frequency				Remarks
		Every 6 Months	Every Year	Every 2 Years	Every 3 Years	
APCI probe						 P.136
Corona needle (for both APCI and DUIS)	-	○				
Interface						 P.159
Spray unit	-	○				
Sampling cone	-	○				Cleaning
Orifice	-			○		Cleaning
Lens system						 P.142
Qarray	-	○				Cleaning
Skimmer	-	○				Cleaning
Multipole 1	-	○				Cleaning
Multipole 2	-	○				Cleaning
Collision cell						 P.30
CID cell ASSY	-		○**			Cleaning



NOTE

Work indicated with ** in the Replacement/Inspection Frequency column must be performed by field engineers.

The cleaning frequency of the interface, lens system, or collision cell may differ depending on the operating conditions at your site.




Fig. 7-1

7.2 Maintenance of the ESI Probe

7.2.1 Replacing the Tubing

A high voltage (up to 5 kV) is applied to the sample solution sprayed from the tip of the capillary ASSY in order to ionize the sample. Since a very fine-gauge tube is used in the capillary ASSY, if any very small pieces of dirt get into the mobile phase or the sample, they can cause a blockage which will lead to a rise in the pump pressure.

If this happens, replace the tubing of the capillary ASSY.

Parts used

Part Name	Part No.
Capillary ASSY	S225-14948-91
	
O-ring	S036-19004-05



NOTE

When carrying out maintenance work on the ESI probe, rest it on the stand.




Fig. 7-2

- 1 Remove the ESI probe from the instrument.



Reference

["3.5.2 Removing the ESI Probe" P.58](#)

- 2 Remove the ESI coupling.
 - 1 Loosen the male nut ① by hand and remove it.
 - 2 Loosen the ESI coupling ②.

3 Remove the knurled screw ③.




Fig. 7-3

4 Remove the ESI coupling ④.





Fig. 7-4

3 Replace the capillary ASSY.

- 1 Remove the capillary ASSY ② from the ESI coupling ①.
- 2 Mount the new capillary ASSY ③ on the ESI coupling ①.
Tighten it firmly by hand.



NOTE

Make sure that the tip of the ESI pipe is not either projecting too far or too recessed.




Fig. 7-5

4 Mount the ESI coupling in the ESI probe.

- 1 Insert the ESI coupling ① into the ESI probe.



Hint

If the O-ring of the ESI coupling is damaged, replace it.



Reference

["10.3.1 ESI Probe: S225-14949-41" P.232](#)

- 2 Tighten the ESI coupling ① by hand and then turn it with a spanner, and adjust the projection of the tip.
Adjust so that the tip of the capillary ASSY projects about 0.5 to 1 mm.




Fig. 7-6

5 Secure the ESI coupling in place.

- 1 Turn the knurled screw **1** by hand and fix the position of the capillary ASSY.
 - 2 Turn the male nut **2** by hand to secure it.
- When securing the male nut **2**, press against the PEEK tube **3** while tightening it.




Fig. 7-7




Fig. 7-8

7.2.2 Replacing the PEEK Tube

As a result of contamination in the mobile phase or sample, the tube may become blocked, causing the pump pressure to rise.

Parts used

Part Name	Part No.
PEEK tube $\phi 1.6 \times \phi 0.13$, red, 3 m	-

1 Remove the ESI probe from the instrument.



Reference
"3.5.2 Removing the ESI Probe" P.58

2 Replace the PEEK tube.

1 Remove the male nuts ❶ at both ends and remove the PEEK tube ❷.

2 Replace the PEEK tube with the new one.

The standard tube length is 420 mm.

3 Attach the tube to the ESI probe and tighten the male nut at each end by hand.




Fig. 7-9



NOTE

- If you are going to cut the PEEK tube before using it, use a tube cutter.
- When securing the male nuts, press against the PEEK tube while tightening them.



Hint

- The tube cutter (S228-32930-01) is not provided as an accessory. It is an option.
- Make sure that the cut face of the PEEK tube is always square. Cutting the tube on an incline will create dead volume which will lead to broadening of peaks.




Fig. 7-10

7.2.3 Replacing the Tube Joint

If the tube joint for nebulizer gas develops a gas leak, replace it.

Parts used

Part Name	Part No.
Tube joint (Elbow, KQ2L23-M5)	S035-60724-21

- 1 Remove the ESI probe from the instrument.



Reference

"3.5.2 Removing the ESI Probe" P.58

- 2 Replace the tube joint.

- 1 Loosen the tube joint with a spanner and remove it.
- 2 Secure the new tube joint with the spanner.




Fig. 7-11



NOTE

Note that overtightening will damage the thread, so care is required.

Finger-tighten the tube joint and then turn it another quarter of a turn with the spanner.

7.3 Maintenance of the APCI Probe

 WARNING	
	The inside of the APCI probe reaches high temperatures and you must wait until the temperature has decreased sufficiently before starting work.
Instructions	Do NOT touch the metal faces at either end or the ceramic cylinder. You could sustain burns.

 CAUTION	
	When performing APCI analysis in which a halogenated mobile phase additive is used, note the following:
Instructions	<p>In the APCI method, the analyte is heated to high temperatures. If the mobile phase contains a halogenated compound such as chloroform, a corrosive gas will be produced. Even a small percentage concentration of halogens in the mobile phase can cause corrosion.</p> <p>Use SUS tubes to connect the high performance liquid chromatograph to the instrument. PEEK resin tubes do not have sufficient strength and should they rupture, solvent may blow out.</p> <p>Black particles may build up on the APCI probe or inside the probe holder.</p> <p>Because of exposure to corrosive gases, you may need to replace the following parts at a shorter interval than specified.</p>

Part Name	Part No.
APCI pipe ASSY	S225-15845-91
Nebulizer joint ASSY	S225-15788-91
Heater unit ASSY	S225-15619-41
Adapter	S225-04993
Ferrule (LCMS-2020)	S225-03748-03
Nut	S225-15739
Heater flange	S225-15486-91
Needle unit	S225-15921-92
DL	S225-15718-91

7.3.1 Replacing the APCI Pipe

Parts used

Part Name	Part No.
APCI pipe ASSY	S225-15845-91
	



NOTE

When carrying out maintenance work on the APCI probe, rest it on the stand.




Fig. 7-12

■ Removing the APCI pipe

- 1 Remove the APCI probe from the instrument.



Reference

["3.6.3 Removing the APCI Probe" P.66](#)

- 2 Remove the APCI cover.

- 1 Remove the screws **1** at three locations.
- 2 Loosen the locking screw **2**.

3 Remove the APCI cover ③.




Fig. 7-13

3 Remove the APCI pipe ASSY.

- 1 Hold the nut ① at the heater side steady and turn the APCI adapter nut ②.
- 2 Remove the APCI pipe ASSY ③.




Fig. 7-14

4 Remove the APCI adapter from the APCI pipe ASSY.




Fig. 7-15



NOTE

Handle the pipe section of the APCI pipe ASSY with care since it is easily bent.

■ Mounting the APCI pipe

1 Mount the APCI pipe.

1 Mount the APCI adapter ② on the new APCI pipe ①.

2 Insert the APCI pipe into the APCI probe.

Insert the tip of the APCI pipe while keeping the pipe vertical.




Fig. 7-16



NOTE

Finger-tighten the APCI adapter and then turn it another quarter of a turn with the spanner.

2 Mount the APCI cover.

1 Tighten the APCI pipe nut ①.

2 Pass the APCI adapter through the hole in the APCI cover ②.

3 Orient the notch ③ correctly.

4 Insert the three projections ④ into the three slots ⑤.




Fig. 7-17

- 5 Secure the APCI cover with the screws ⑥ at three locations.
- 6 Tighten the locking screw ⑦.




Fig. 7-18

7.3.2 Replacing the Heater Unit

Parts used

Part Name	Part No.
Heater unit ASSY	S225-15619-41

■ Removing the APCI heater

- 1 Remove the APCI probe from the instrument.
 - Reference
["3.6.3 Removing the APCI Probe" P.66](#)
- 2 Remove the APCI heater from the APCI probe.
 - 1 Remove the screws ① at three locations.
 - 2 Loosen the locking screw ②.
Use a hexagon wrench 1.5.




Fig. 7-19

- 3 Remove the APCI cover ③.
- 4 Remove the screws ④ at three locations.
- 5 Draw out the heater unit ⑤.




Fig. 7-20

3 Remove the nebulizer unit.

- 1 Remove the screws ① at two locations.
 - 2 Remove the nut.
- Hold the nut at the nebulizer side steady and loosen the nut at the heater side ②.




Fig. 7-21

3 Remove the nebulizer unit.




Fig. 7-22

■ Mounting the APCI heater

1 Mount the heater in the nebulizer unit.

1 Mount the ferrule **1**, nut **2** and heater **3**.

2 Tighten the two screws **4**.

Tighten these screws in order to achieve alignment with the position of the ferrule.




Fig. 7-23



NOTE

Take care to mount the ferrule in the correct orientation.

2 Secure the nut.

Hold the nut at the nebulizer side steady while tightening the nut at the heater side.




Fig. 7-24



NOTE

Push the heater toward the nut when tightening the nut at the heater side.

3 Insert the heater unit.

- 1 Align the cable section of the heater unit with the recess ❶ in the probe pedestal.
- 2 Align the tip of the heater unit with the hole in the probe pedestal.




Fig. 7-25

4 Mount the APCI cover.

- 1 Tighten the screws ❶ at three locations.
- 2 Pass the adapter through the hole in the APCI cover ❷.
- 3 Orient the notch ❸ correctly.
- 4 Insert the three projections ❹ into the three slots.




Fig. 7-26

- 5 Tighten the screws ❺ at three locations and the locking screw ❻.




Fig. 7-27

7.3.3 Cleaning and Replacement of the APCI Corona Needle

WARNING



Instructions

Before starting maintenance work, turn the high-voltage switch OFF in the LabSolutions program and disconnect the high-voltage cable.

If the high-voltage cable is not disconnected there will be a danger of electric shock.



CAUTION



Instructions

Take care because the tip of the APCI needle is sharp.

Parts used


Part Name	Part No.
APCI needle	S225-15877-92
	
Needle unit ASSY	S225-14290-41
	

1 Remove the needle unit from the instrument.



Reference

["3.6.4 Removing the Corona Needle for APCI" P.69](#)



If the tip of the needle is blackened or there is a foreign body present, clean the needle.

Fig. 7-28

2 Polish the APCI needle.

1 Remove the APCI needle ❶ with a screwdriver.

2 Polish the tip ❷ of the APCI needle.

Polish the tip of the APCI needle using about 4-micron lapping film. When polishing, take care not to curve the tip.




Fig. 7-29

3 Immerse the APCI needle in methanol and clean it by ultrasonic cleaning.



Fig. 7-30

4 Mount the APCI needle on the support arm.




Fig. 7-31

5 Mount the needle unit on the instrument.

Reference

["3.6.1 Mounting the Corona Needle for APCI" P.60](#)

NOTE

Before mounting the needle unit, check that the tip of the APCI needle is not bent by using the needle alignment tool.

Check that the needle is contained within the range of the slit and repair the needle if it is bent.




Fig. 7-32

7.3.4 Replacing the Gas Joint

If the tube joint for the nebulizer gas develops a gas leak it must be replaced.

Parts used

Part Name	Part No.
Gas joint (Half union, KQ2H01-M5)	S035-60725-01

1 Remove the APCI probe from the instrument.

Reference

["3.6.3 Removing the APCI Probe" P.66](#)

- 2 Loosen the tube joint with a spanner and remove it.




Fig. 7-33

- 3 Secure the new tube joint with the spanner.




Fig. 7-34

**NOTE**

Finger-tighten the tube joint and then turn it another quarter of a turn with the spanner. Note that overtightening will damage the thread, so care is required.



7.4 Maintenance of DUIS

For maintenance procedures relating to the ESI probe, see "3.7 Preparation for DUIS Analysis" P.72.

 CAUTION	
	Take care because the tip of the DUIS needle is sharp.
Instructions	

7.4.1 Cleaning or Replacing the Needle for DUIS

Parts used


Part Name	Part No.
DUIS needle	S225-15877-92
	
Needle unit ASSY	S225-14290-41
	

- 1 Remove the needle unit from the instrument.



Reference

["3.7.3 Removing the ESI Probe" P.76](#)



If the tip of the needle is blackened or there is a foreign body present, clean the needle.

Fig. 7-35

2 Polish the DUIS needle.

- 1 Loosen the DUIS needle ❶ with a screwdriver and remove it.
- 2 Polish the tip ❷ of the DUIS needle.

Polish the tip of the needle using about 4-micron lapping film. When polishing, take care not to curve the tip.




Fig. 7-36

3 Immerse the DUIS needle in methanol and clean it by ultrasonic cleaning.



Fig. 7-37

4 Mount the DUIS needle on the support arm.




Fig. 7-38

NOTE

Check the mounting angle of the DUIS needle with the tool and, if the angle is incorrect, correct it.

Reference

["7.3.3 Cleaning and Replacement of the APCI Corona Needle" P.136](#)

5 Mount the needle unit to the instrument.

Reference

["3.7.1 Mounting the Corona Needle for DUIS" P.72](#)

7.5 Maintenance of the Lens System

WARNING



Stop the vacuum system and turn the instrument power switch OFF before starting maintenance work.

Instructions

If you do not turn the power switch OFF there will be a danger of electric shock.

CAUTION



Wear clean gloves when carrying out maintenance work on the lens system.

Instructions

Since the lens system focuses the ions generated by the ionization unit, it becomes soiled over long periods of use and the detection sensitivity diminishes. For this reason you must clean the following lens system components: Qarray, skimmer, multipoles, and entrance lens. Also, clean the lens system components when a considerable decrease in sensitivity is observed.

7.5.1 Removing the Multipoles and Entrance Lens

- 1 Stop evacuation and turn the power switch OFF.



Reference

["3.2 Stopping the Instrument" P.48](#)

- 2 Open the multipole maintenance cover.
 - 1 Press the indentation with your finger.
 - 2 Lift the entire cover upward to remove it.




Fig. 7-39

3 Open the lens system door.

Loosen the knurled screw to open the door. The vacuum housing must be opened to atmosphere before opening the door.




Fig. 7-40 (The figure shows the LCMS-8040.)

4 Remove multipole 2 and the entrance lens.

- 1 Remove the hooks at two locations from the hook plate.
- 2 Slide the flange section of multipole 2 and lift it upward.




Fig. 7-41 (The figure shows the LCMS-8040.)

Hint

Once you have removed the flange section, lift the multipole upward while keeping it horizontal.

NOTE

- Hold the flange section during removal so that you do not touch the electrode unit.
- Do not place any parts on top of the instrument.

- 5 Remove multipole 1 and the entrance lens.
- 1 Remove the hooks at two locations from the hook plate.
 - 2 Slide the flange section of multipole 1 and lift it upward.




Fig. 7-42 (The figure shows the LCMS-8040.)

7.5.2 Cleaning the Multipoles and Entrance Lens

CAUTION	
	Do NOT disassemble either multipole body. Disassembling them could cause injuries and instrument failure.
Prohibitions	

- 1 Clean the multipoles.
 - 1 Immerse the multipoles in methanol and clean them by ultrasonic cleaning.
 - 2 Remove the multipoles from the methanol and leave them to air until dry.




Fig. 7-43

**NOTE**

Do not use paper towels but allow the multipoles to dry naturally.

7.5.3 Removing the Qarray and Skimmer

- 1 Open the probe cover.
- 2 Open the source window and remove it.




Fig. 7-44

- 3 Remove the ESI and APCI probes from the instrument.



Reference

["3.5.2 Removing the ESI Probe" P.58](#)

["3.6.3 Removing the APCI Probe" P.66](#)

- 4 Open the front door.
- 5 Pull the lever up and open the probe holder.

**NOTE**

The vacuum housing must be opened to atmosphere before opening the door.




Fig. 7-45

6 Remove the flange.




Fig. 7-46

7 Draw out the Qarray and skimmer.




Fig. 7-47

7.5.4 Cleaning the Qarray and Skimmer

CAUTION



Prohibitions

Do NOT disassemble the Qarray body.
Disassembling it could cause injuries and instrument failure.



NOTE

Rest the Qarray on the probe stand and then remove the skimmer.




Fig. 7-48

1 Clean the Qarray and skimmer.

- 1 Remove the two screws **1** and remove the skimmer **2**.
- 2 Immerse the skimmer and Qarray in methanol and subject them to ultrasonic cleaning.




Fig. 7-49

2 Return the parts to their original positions.

1 Dry the parts after cleaning.



NOTE

If the central part of the inner face of the skimmer is soiled, wipe it clean, e.g. with a cloth. If there is any contamination adhering to the central hole, remove it.




Fig. 7-50

2 Mount the skimmer ❶ and secure it with the two screws ❷.




Fig. 7-51

7.5.5 Mounting the Qarray and Skimmer

- 1 Mount the Qarray and skimmer by fitting them on the guide pins.
Insert the Qarray all the way in until it slightly pushes back towards you on the springs.




Fig. 7-52

- 2 Attach the flange.




Fig. 7-53

- 3 Return the probe holder to its original position and close the lever.




Fig. 7-54

- 4 Mount the ESI and APCI probes.

 Reference

["3.7.2 Mounting the ESI Probe" P.75](#)

["3.6.2 Mounting the APCI Probe" P.63](#)

- 5 Mount the source window.

7.5.6 Mounting the Multipoles and Entrance Lens

1 Mount the multipoles and the entrance lens.

- 1 Hold and slide the flange section and insert multipole 1.
- 2 Attach the hooks at two locations onto the hook plate.




Fig. 7-55 (The figure shows the LCMS-8040.)

2 Mount multipole 2 and the entrance lens.

- 1 Hold the flange section and insert multipole 2.
- 2 Attach the hooks at two locations onto the hook plate.




Fig. 7-56 (The figure shows the LCMS-8040.)

- 3 Close the lens system door.
Close the lens system door and secure it with the knurled screw.




Fig. 7-57

7.5.7 Replacing the O-Ring

- 1 Stop evacuation and turn the power switch OFF.



Reference

["3.2 Stopping the Instrument" P.48](#)

- 2 Open the multipole maintenance cover.
 - 1 Press the indentation with your finger.
 - 2 Lift the entire cover upward to remove it.

- 3 Loosen the knurled screw and open the lens system door.




Fig. 7-58 (The figure shows the LCMS-8040.)

- 4 Remove the O-ring and, if there is any damage, replace it.




Fig. 7-59 (The figure shows the LCMS-8040.)



NOTE

When inserting the O-ring into the groove, make sure that there is no dirt adhering to it.

7.6 Maintenance of the Standard Sample Introduction Unit

7.6.1 Replacing the Standard Sample

- 1 Stop pumping the standard sample.

 Reference

"8.2.1 Checking Ions in the [MS Tuning] Window" P.201

- 2 Open the front door.




Fig. 7-60

- 3 Remove the standard sample bottle.

- 1 Disconnect the resistance tube.

- 2 Lift the standard sample bottle up to remove it from the bracket.

- 3 Turn the bottle ② while holding the bottle cap ①.




Fig. 7-61

 **NOTE**

When removing the standard sample bottle, take care not to injure yourself on the edge of the bracket.

 **Hint**

Some packing may become stuck in the bottle cap.

4 Pour standard sample into the bottle.




Fig. 7-62

NOTE

If there is only a little standard sample left in the bottle, throw it away and wash the bottle before pouring in new standard sample. Fill the bottle with at least 40 to 80 mL of standard sample (the highest reading on the bottle's scale corresponds to 80 mL).

Reference

["8.3 Standard Sample" P.213](#)

5 Return the standard sample bottle to its original position.

- 1 Attach the bottle cap **1** to the bottle **2**.
- 2 Fit the standard sample bottle onto the bracket.
- 3 Attach the resistance tube **3**.




Fig. 7-63

NOTE

- Always ensure there is packing between the bottle and the bottle cap.
- When mounting the standard sample bottle, take care not to injure yourself on the edge of the bracket.

6 Close the front door.

7.6.2 Replacing the Resistance Tube

- 1 Stop pumping the standard sample.



Reference

"8.2.1 Checking Ions in the [MS Tuning] Window" P.201

- 2 Open the front door.




Fig. 7-64

- 3 Replace the resistance tube.
 - 1 Disconnect the resistance tube ❶.
 - 2 Remove the male nuts ❷ from the resistance tube ❶.
 - 3 Fit the male nuts ❷ to the new resistance tube ❸.




Fig. 7-65

- 4 Fit the resistance tube between the standard sample bottle and the connector.



NOTE

Take care not to bend the resistance tube too much.

Maintain a moderate bending radius of at least 40 mm while handling the tube.

- 5 Close the front door.

7.6.3 Replacing the Filter

- 1 Remove the standard sample bottle from the instrument.



Reference

["7.6.1 Replacing the Standard Sample" P.154](#)

- 2 Remove the gas tubing for delivering the standard sample from the standard sample bottle.

The tube can be removed by pressing the release bush of the gas one-touch joint.




Fig. 7-66

- 3 Remove the two screws and disassemble the parts.




Fig. 7-67

4 Secure the base with a spanner.




Fig. 7-68

5 Replace the filter with a new one and tighten the plug.




Fig. 7-69

NOTE

To tighten the plug, first hand-tighten it and then turn it another 60 to 90 degrees.

6 Tighten the two screws to assemble the parts.




Fig. 7-70

7.7 Cleaning the Spray Unit

WARNING



Instructions

- Before starting maintenance work, turn the heater OFF from the LabSolutions program and make sure that the temperature of the heated block has fallen to 50 °C or lower.
The spray unit reaches high temperatures and could cause burns.
- Before starting maintenance work, turn the high-voltage switch OFF in the LabSolutions program and disconnect the high-voltage cable.
If the high-voltage cable is not disconnected there will be a danger of electric shock.

If soiling cannot be removed, replace the part.

Parts used

Part Name	Part No.
Sampling cone	S225-15487
Flat countersunk head screw M3 × 5	S225-14287-41
Drain mesh	-

1 Remove the ionization probe.

- 1 Open the probe cover.
- 2 Open the source window.
- 3 Remove the ionization probe.

With APCI or DUIS, remove the corona needle.




Fig. 7-71

Reference

["3.5.2 Removing the ESI Probe" P.58](#)

["3.6.3 Removing the APCI Probe" P.66](#)

["3.6.4 Removing the Corona Needle for APCI" P.69](#)

■ Routine inspection

1 Wipe off soiling.

Moisten gauze with a solvent that can remove soiling (such as water/methanol) and wipe the soiling off with the gauze.

Clean the locations indicated below.

- Sampling cone
- Heated block
- Heater flange
- Inside face of the source window
- Inside walls of the spray unit




Fig. 7-72

NOTE

Do not clean the tip of the DL since this can cause blockage of the tubing. If the tip of the DL is badly soiled, replace the DL.

Reference

["7.9.1 Removing the DL from the Instrument" P.163](#)

■ If the soiling is severe:

1 Clean the parts.

- 1 Remove the two flat countersunk head screws ❶ and remove the sampling cone ❷.



NOTE

Remove the flat countersunk head screws only after the temperature drops to around room temperature in order to prevent damage to the screws.

- 2 Remove the drain mesh ❸.

Lift it up with tweezers.

- 3 Wipe the rear of the sampling cone.

Moisten gauze with a solvent that can remove soiling (such as water/methanol) and wipe the soiling off with the gauze.

If the heated block is soiled, scrub it with a nylon brush and then wipe it clean with gauze moistened with a solvent that can remove the soiling (such as water/methanol).




Fig. 7-73

- 4 Subject the sampling cone, flat countersunk head screws and drain mesh to ultrasonic cleaning in a methanol solution.



HINT

If the soiling is particularly bad, polish the parts with about 4-micron lapping film.




Fig. 7-74

- 2 After drying the parts, mount them in the instrument.



NOTE

When mounting the parts, check if there is any foreign body in the hole of the sampling cone, and if there is any cleaning fluid left on the parts. If there is a foreign body it may cause blockage of the DL.

7.8 Replacing the O-Ring of the Source Window

If the O-ring of the source window has become damaged or has deteriorated, it could lead to solvent leakage.

In such cases it must be replaced.

Parts used

Part Name	Part No.
O-ring, 4D G145	S036-12527

- 1 Remove the O-ring.
 - 1 Open the probe cover ①.
 - 2 Open the source window ②.
 - 3 Remove the O-ring ③.




Fig. 7-75



NOTE

If the O-ring is difficult to remove, remove the source window and then remove the O-ring. The source window can be removed by lifting it upward.

- 2 Fit the new O-ring.



NOTE

When fitting the O-ring, be careful not to damage it or let any dirt adhere to it.

7.9 Replacing the DL

WARNING



Instructions

- Before starting maintenance work, turn the heater OFF from the LabSolutions program and make sure that the temperature of the heated block has fallen to 50 °C or lower.
The spray unit reaches high temperatures and could cause burns.
- Before starting maintenance work, turn the high-voltage switch OFF in the LabSolutions program and disconnect the high-voltage cable.
If the high-voltage cable is not disconnected there will be a danger of electric shock.

Parts used

Part Name	Part No.
DL ASSY	S225-15718-91

7.9.1 Removing the DL from the Instrument

- 1 Open the probe cover **1**.
- 2 Unlock and open the source window **2**.
Grip the lower right of the source window and pull towards you.
- 3 Remove the ionization probe **3**.
If using APCI or DUIS, remove the corona needle too.
- 4 Loosen the two screws **4** with the hexagon wrench provided as an accessory.




Fig. 7-76

- 5 Insert the drawing tool **5** provided as an accessory under the heated block and remove the heater flange **6**.




Fig. 7-77

- 6 Disconnect the connector **7**.
Press in the catches on the connector as you remove it.
- 7 Loosen the two DL locking screws **8** with the hexagon wrench provided as an accessory.
Turn the hexagon wrench about three turns.
- 8 Turn the notched section **9** to the direction in which the DL is to be removed, then draw it out.




Fig. 7-78

**Hint**

The triple inlet turbo molecular pump changes to standby mode when the DL is removed. The pump returns to normal mode about 3 minutes after mounting the DL (connector **7**).

7.9.2 Mounting the DL in the Instrument

- 1 Turn the notched section in the direction in which it engages to secure the DL to the heater flange ④, then connect the connector ①.
- 2 Tighten the two DL locking screws ② with the hexagon wrench provided as an accessory.
Tighten the locking screws sufficiently.
- 3 Insert the heater flange ④ into the IF flange and tighten the two screws ⑤.
Use the drawing jig to insert the heater flange.




Fig. 7-79



NOTE

When mounting the heater flange, take care about the positioning of the top and bottom of the flange.

- 4 Mount the ionization probe ①.
- 5 Close the source window ②.
- 6 Close the probe cover ③.




Fig. 7-80

7.10 Cleaning and Replacing the Orifice

 WARNING	
 Instructions	<ul style="list-style-type: none"> • Stop the vacuum system and turn the instrument power switch OFF before starting maintenance work. <p>If you do not turn the power switch OFF there will be a danger of electric shock.</p> <ul style="list-style-type: none"> • Before starting maintenance work, turn the heater OFF from the LabSolutions program and make sure that the temperature of the heated block has fallen to 50 °C or lower. <p>The spray unit reaches high temperatures and could cause burns.</p> <ul style="list-style-type: none"> • Before starting maintenance work, turn the high-voltage switch OFF in the LabSolutions program and disconnect the high-voltage cable. <p>If the high-voltage cable is not disconnected there will be a danger of electric shock.</p>

Parts used

Part Name	Part No.
Orifice	S225-15479

7.10.1 Removing the Orifice from the Instrument

- 1
 Stop the vacuum.
 -  Reference
["3.2.1 Stopping the Vacuum System" P.48](#)
- 2
 Remove the DL.
 -  Reference
["7.9.1 Removing the DL from the Instrument" P.163](#)

3 Remove the orifice ASSY.

- 1 Loosen the single screw ❶.
- 2 Turn the notched section ❷ in the direction indicated by the arrow and remove the orifice ASSY.




Fig. 7-81

4 Turn the orifice with a spanner and remove it.




Fig. 7-82

7.10.2 Cleaning and Assembling the Orifice

- 1 Use approximately 4-micron lapping film to clean and polish the tip of the orifice until you can see its metal surface.




Fig. 7-83

- 2 Wipe off any soiling on the orifice tip using solvent. Moisten gauze with a solvent that can remove soiling (such as water/methanol) and wipe any soiling off with the gauze.




Fig. 7-84

- 3 Dry the orifice parts.

4 Check that there are no foreign bodies in or around the orifice.




Fig. 7-85

**NOTE**

If there is a foreign body in the orifice, use the provided "DL plug" and the air duster to remove it.

5 Turn the orifice with a spanner to tighten it.




Fig. 7-86

**NOTE**

Refit the O-ring if it has come loose from the orifice. Replace the O-ring if it is damaged in any way.

7.10.3 Mounting the Orifice in the Instrument

- 1 Mount the orifice ASSY.
 - 1 Turn the notched section ❶ in the direction indicated by the arrow to mount the orifice ASSY.
 - 2 Tighten the single screw ❷.




Fig. 7-87

**NOTE**

If the O-ring (P/N: S036-19004-11) is damaged, replace it.

- 2 Mount the DL.
 -  Reference
["7.9 Replacing the DL" P.163](#)
- 3 Start the vacuum.
 -  Reference
["3.1.2 Starting the Vacuum System" P.42](#)

7.11 Maintenance of the Rotary Pump



WARNING



Instructions

Stop the vacuum system and turn the instrument power switch OFF before starting maintenance work.

If you do not turn the power switch OFF there will be a danger of electric shock.

■ Rotary pump oil change

The rotary pump (type E2M28) that is used as the auxiliary pump for this product's vacuum system requires an oil change every four months. Failure to change the oil will cause trouble including vacuum deficiency, oil leaks and increased noise; it is therefore essential to perform oil changes.

Use an oil whose characteristics are matched to the model of the pump and which is authorized by Shimadzu.

Parts used

Part Name	Part No.
Rotary pump oil Ultragrade19 (4 L container)	S017-30163-02



WARNING



Instructions

Change the oil in response to change in its color and in its volume.

It may be necessary to change the oil more frequently than at four-month intervals, depending on the conditions of use and frequency of use of the instrument.

7.11.1 Changing the Oil



NOTE

Immediately after stopping the instrument the rotary pump oil will be at a high temperature (approximately 30 °C above room temperature), so you must wait 10 minutes before changing the oil.

- 1 Stop the instrument.



Reference

["3.2 Stopping the Instrument" P.48](#)

- 2 Wait about 10 minutes after stopping it.

- 3 Inject the new oil.

- 1 Turn the motor switch of the rotary pump OFF.




Fig. 7-88

- 2 Remove the drain plug.

Drain the oil while receiving it in a tray or plastic bag with a capacity of about 2 L. Note that oil may spray out when the drain plug is removed, so care is required.

- 3 When oil has stopped coming out of the drain port, close the drain plug.

- 4 Remove the oil filler plug.

- 5 Pour in new oil up to the bottom of the MAX indication of the oil gauge (approximately 1.5 L).





Fig. 7-89

- 4 Close the oil filler plug.
- 5 Turn the motor switch of the rotary pump ON.
- 6 Reset the rotary pump oil change frequency.

- 1 Click the  (System Check) icon.




Hint


Alternatively, select [Instrument] - [System Check...].

The [System Check] window is displayed.

- 2 Click [Reset...].
If [Reset...] is not displayed, click [Advanced >>].



- 3 Enter the rotary pump oil use time (set to 0 after an oil change) and click [Set].



- 4 Click [Close].

7.11.2 Gas Ballast Valve

CAUTION



Instructions

Under normal conditions of use, open the gas ballast valve once a week for around 15 to 30 minutes to remove the contaminants in the oil.

If the water content ratio of the mobile phase is high or the flow volume is large, open the gas ballast valve once every 1 to 3 days.

When you open the valve, a large quantity of oil mist will emerge from the exhaust port. Make sure that the exhaust gas is always released into a duct system.

Since a large quantity of mobile phase and water is taken into the rotary pump in the instrument, using it over a long period of time without opening the gas ballast valve will cause deterioration in exhaust performance and instrument failure.

1 Exit LabSolutions.



Reference

"3.2.2 Exiting LabSolutions" P.48

2 Open the gas ballast valve and wait 15 to 30 minutes.




Fig. 7-90

3 Close the gas ballast valve.



NOTE

Apart from the oil changes described above, the rotary pump is also subject to some checks in periodic inspections.

For details on these checks, refer to 5 "Maintenance" in the separately attached E2M28 instruction manual (A373-10-880).

7.12 Vacuum Leak Check

Running the instrument while there is a vacuum leak will lead to problems such as reduced sensitivity, increased noise, and breakage of the ion gauge (vacuum gauge) filament.

For this reason, always carry out a vacuum leak check on restarting the instrument after it has been stopped.

A whole range of vacuum leaks can affect the MS, from large leaks that stop the rotary pump from functioning normally to small leaks that are no hindrance to the normal operation of the triple inlet turbo molecular pump. These vacuum leaks develop from very small problems.

There are always O-rings fitted on parts with connections, such as the probe holder, between the lens system door and its housing, as well as the adjacent parts of the orifice. By making close contact with both sides at connections, these O-rings maintain the vacuum inside the housing by preventing vacuum leaks. If dirt adheres to an O-ring, or if you forget to fit an O-ring, a vacuum leak will occur.

When you stop the instrument and carry out maintenance work with the vacuum applied (cleaning of the lens system, replacement of the DL, replacement of the orifice), observe the following points in order to prevent vacuum leaks.

Cautions for preventing vacuum leaks

 CAUTION	
 Prohibitions	<ul style="list-style-type: none"> • Do NOT damage the seal faces. • Do NOT forget to fit O-rings.
 Instructions	<ul style="list-style-type: none"> • Fit O-rings correctly. • Remove any dirt that adheres to O-rings. • Remove dirt on the faces with which the O-ring makes close contact (seal faces).

In the following circumstances, the rotary pump and triple inlet turbo molecular pump are not operating normally.

- The exhaust noise of the rotary pump continues for several minutes.
- The triple inlet turbo molecular pump started, but after the elapse of around 10 minutes the power to the vacuum system went OFF automatically.

In these cases it is probable that there is a vacuum leak or that the vacuum pump is abnormal. After the vacuum housing has leaked to atmospheric pressure, check the probe holder, lens system door, DL and orifice in accordance with "[Cautions for preventing vacuum leaks](#)" described previously.

If it is thought that a vacuum leak has occurred for reasons other than those described above, contact your Shimadzu representative.

NOTE

If the pressure is high, focus your attention on the locations subjected to maintenance.

Parts where there is a possibility of vacuum leakage




Fig. 7-91

7.13 Repairing the Nebulizer Gas Tube

- 1 Open the probe cover.
- 2 Press the release bush and remove the nebulizer gas tube.
Stop the gas supply from LabSolutions before starting work.




Fig. 7-92

- 3 Cut perpendicularly across the tube about 1 cm from the end with a cutter knife.
Alternatively, cut it with the tube cutter (option P/N S228-32930-01).




Fig. 7-93

- 4 Fit the nebulizer gas tube by inserting it as far as it will go.

7.14 Replacing the Waste Tubes

Waste tubes are fitted at the locations indicated below.

If these tubes are particularly soiled, or are damaged, replace them.




Fig. 7-94



NOTE

Replace the waste tubes once every three years.

7.15 Cleaning the Leak Tray

The leak tray is fitted with a leak sensor to detect liquid leakage. If a liquid leakage occurs, wipe up all the liquid in the leak tray by following the procedure below.

- 1 Open the front door.




Fig. 7-95

- 2 Completely wipe up all leakage around the leak sensor.
If there is any soiling, wipe it off with water.




Fig. 7-96

7.16 Cleaning the Exterior

If the body cover is soiled, wipe it with a soft dry cloth or tissue paper.

If the soiling is particularly bad, clean it off as follows.

- 1 Wipe the cover with a cloth that has been dipped in a diluted neutral detergent and wrung out well.
- 2 Dip a cloth in water, wring it out well and wipe the body so that no detergent remains, then wipe off the moisture with a dry cloth.



NOTE

Do not leave the body moist with water or wipe it over using any kind of alcohol or thinner. Doing so could cause rust formation or discoloration.

7.16.1 Cleaning the Fans

- 1 Remove the dust from the two fan filters on the instrument with, for example, a vacuum cleaner. Clean the fan filters about every six months.




Fig. 7-97

This page is intentionally left blank.

8.1 Installation

8.1.1 Conditions for the Installation Site

The instrument will be installed and adjusted by your Shimadzu representative.

In order to keep the instrument operating stably and to obtain highly reliable analysis results, observe the conditions given below.

■ Power supply

 WARNING	
 Prohibitions	<ul style="list-style-type: none">• Do NOT share the power supply for this instrument with other devices.• The electric capacity of the instrument must not exceed the stipulated capacity of the earth leakage breaker.• Do NOT rest anything heavy on the power cable. Do NOT run the power cable close to any heat-generating equipment. Such treatment could damage the cable, leading to fire, electric shock or instrument failure. If the cable becomes damaged, contact your Shimadzu representative immediately.
 Instructions	<ul style="list-style-type: none">• Connect the instrument to a compatible power supply. This product has a power supply voltage of 230 V AC and a power consumption of 3.45 kVA. Connect to a compatible power supply. A supply that is not compatible will cause fire or electric shocks. If the power supply voltage is unstable or the power supply capacity is insufficient, the intended performance of the instrument cannot be guaranteed. You must check the power supply capacity of the system as a whole and prepare an appropriate supply.• The power supply that you connect to must be equipped with an earth leakage breaker.• The instrument must be grounded (earthed). If you do not ground it, there will be a risk of instrument failure and electric shock in the event of current leaks. It is also essential to ground the instrument in order to stabilize it.

(1) Required power supply

LCMS-8030/LCMS-8040	Single phase 230 V AC, 15 A (50/60 Hz)
Power supply voltage range to guarantee intended performance:	218.5 to 241.5 V
Power supply voltage range to guarantee operation:	207 to 253 V
Frequency stability:	± 0.5 Hz max.
(200 V max.)	With APCI 10 A (200 V, 50/60 Hz) Without APCI 9 A (200 V, 50/60 Hz)
(200 V steady)	4 to 6 A (200 V, 50/60 Hz)
(230 V max.)	With APCI 10 A (230 V, 50/60 Hz) Without APCI 9 A (230 V, 50/60 Hz)
(230 V steady)	4 to 5 A (230 V, 50/60 Hz)
LC	Refer to the specifications or the instruction manual provided with your LC.

(The electric capacity varies depending on the configuration of the LC system.)

Prepare the power supply for the LC in accordance with the specifications of the LC unit you are using or the instruction manual.

For details on the power supply specifications of peripheral devices such as the computer and printer, check the specifications or instruction manuals of the units you are using. Approximate values are indicated below.

Power for PC system	Supply correct voltage to the PC or printer, etc. referring to the instruction manuals accompanying the PC or printer, etc.
---------------------	---

**CAUTION**

Use a PC, printer, etc. that conform to your local laws and regulations.

Instructions

(2) Grounding

100 Ω max.

In order to prevent electric shocks, be sure to ground the instrument (connect it to earth).



(3) Installation conditions

OVERVOLTAGE CATEGORIE: CATEGORIE II (IEC)

POLLUTION DEGREE: 2 (IEC)

■ Installation environment

 WARNING	
 Prohibitions	<p>Do NOT use naked flames.</p> <p>The use of naked flames is prohibited at the site where the high performance liquid chromatograph mass spectrometer is installed. You should also avoid installing other equipment that produces a naked flame in the same room. To ensure readiness in the event of an accident, install a fire extinguisher.</p>
 Instructions	<ul style="list-style-type: none"> • Be sure to provide ventilation in the room. Some of the solvents that are used with a high performance liquid chromatograph mass spectrometer are inflammable or toxic. Also note that this instrument contains a large quantity of nitrogen gas. Its use in a room that is inadequately ventilated could cause oxygen deficiency. <p>Install the instrument in a room that has a ventilation mechanism such as the type of draft chamber in general use (approx. 20 m³/min), and feed the exhaust tube into the draft chamber.</p> • Use a duct system for exhaust. <p>Be sure to release the exhaust gas from the rotary pump, the solvent vapor that builds up in the waste container, and nitrogen gas into a duct system such as a draft chamber. Be sure to provide separate exhaust channels for the exhaust from the rotary pump and the nitrogen for ionization.</p> <p>Failure to do so will lead to contamination of the mass spectrometer.</p> • Install a wash basin. <p>If a solvent gets into someone's eyes or someone touches a toxic solvent, it has to be rinsed away immediately. Install a wash basin as close to this product as possible.</p>



 CAUTION	
 Prohibitions	<ul style="list-style-type: none"> Do NOT install this product in a location where there is corrosive gas, contaminants or a lot of dust. <p>Avoid installing the instrument at a location where there is corrosive gas or a lot of dust in order to maintain its performance.</p> <ul style="list-style-type: none"> Do NOT install the instrument close to any device that generates a strong magnetic field. <p>To ensure that the instrument can be operated normally, do not install it at a location where there is a strong magnetic field.</p> <p>In addition, if there is a lot of noise in the power line add a noise filter.</p>
 Instructions	<p>Observe the installation conditions.</p> <ul style="list-style-type: none"> A room with a temperature within the range 18 to 28 °C and where the change in room temperature throughout the day is small A location where the instrument is not directly exposed to the airflow from a heater/air conditioner A location not exposed to direct sunlight A location where there is little vibration A location where the humidity remains within the range 40 to 70 %

(1) Temperature, humidity

Temperature 18 °C to 28 °C

Humidity 40 % to 70 % (to be no dew condensation)

Note, however, that temperature changes over a short time can cause instability in the instrument. In order to ensure that stable analysis is performed, make sure that changes in temperature over a short period of time do not exceed 3 °C.

(2) Installation space

The table on which the instrument is installed must be robust, steady and flat, and must be able to comfortably bear the weight of the instrument. The body of the instrument weighs approximately 130 kg (excluding the rotary pump). You must also take into account the weight of the LC system used, and the weight of the computer system.

For an example installation see "[9.5 Example Installation](#)" P.222. Note that, for maintenance purposes, a space of at least 30 cm is required at the right side of the LCMS-8030/LCMS-8040 main body. A space of at least 5 cm is also required at the left side of the instrument. Install the instrument so that it is separated by at least 5 cm from the LC unit.

(3) Ventilation

Some of the solvents used with this instrument are inflammable or toxic. The instrument also uses a large quantity of nitrogen gas, so its use in a room with inadequate ventilation could cause oxygen deficiency. Install equipment that provides adequate ventilation in the room.

(4) Exhaust gas from the rotary pump, solvent vapor

Be sure to release the exhaust gas from the rotary pump and the solvent vapor that builds up in the waste container into a duct system such as a draft chamber.

Be sure to provide separate exhaust lines for the exhaust from the rotary pump and the nitrogen for ionization. Failure to do so will lead to contamination of the mass spectrometer.

(5) Other considerations

Install the instrument in a location free from contaminants, dust, vibration, electromagnetic noise, corrosive gas, interfering magnetic fields, and other hindrances.

In order to maintain the instrument's performance, also give adequate consideration to the following points.

1. Changes in room temperature during use must be small.
2. The airflow from heaters, air conditioners, etc. must not directly strike the instrument.
3. The instrument must not be exposed to direct sunlight.

■ Installation space




CAUTION




Instructions

- Install the instrument on a desk or table that meets the following conditions.
 - Flat and stable / Able to comfortably bear the weight of this product (approx. 130 kg) and the entire computer system / Having a depth of at least 75 cm
Failing to meet these conditions can cause accidents in which the table collapses or units fall off.
 - Allowing the instrument to be placed a sufficient distance from walls
Make the distance between the rear face of the instrument and the wall at least 30 cm and the distance between the left and right faces and other units at least 5 cm. If these conditions are not met it will not be possible for the fans to provide adequate air cooling and there will be a risk that the instrument will overheat and its performance will drop.
 - Providing sufficient space for maintenance
For maintenance purposes a space of at least 30 cm is required to the right of the body of the instrument. Place units that are easy to move, like the computer and printer, on the right side of the instrument body, and place walls and units that are difficult to move at a distance of at least 30 cm from the instrument body.
 - Having an appropriate height
Use a table with a height of about 70 to 80 cm. Using a table with a height outside this range will impair the operating convenience of the instrument, so care is required.



- Install the instrument so that the power switch can be operated easily.

If the power switch is difficult to operate, it will not be possible to switch the power off immediately in an emergency.



Power switch

CAUTION

Instructions


- Install the rotary pump on a stable floor.

Install the body of the instrument and the rotary pump at a distance that allows them to be connected with the tubing.


If the tubing is extended the intended performance of the machine cannot be guaranteed.

An example configuration of a typical system and its installation space are shown in Fig. 8-1.

[Plan View]



[Front Elevation]



Units: mm

Fig. 8-1

8.1.2 Connecting the Instrument

■ Connecting the rotary pump

- 1 Connect the vacuum hose of the rotary pump.
- 2 Connect the rotary pump mist exhaust hose.
Connect the other end of the hose to the exhaust system.




Fig. 8-2

- 3 Connect the rotary pump's power cable to the instrument.
- 4 Turn the rotary pump's power switch ON.




Fig. 8-3

■ Fitting the oil mist filter and oil return tube

The oil mist filter is a part that traps the exhaust mist from the rotary pump.

The oil return tube serves to return the oil that accumulates in the oil mist filter.

In cases where the quantity of mobile phase introduced for mass spectrometer is 0.5 mL/min or greater, a large quantity of water is introduced into the rotary pump and for this reason the gas ballast valve has to be opened and closed frequently in order to remove this water. In cases like this you are recommended to use the oil return kit.

The rotary pump oil return kit is included in the optional parts and should be purchased if required.

Part Name	Part No.
Rotary pump oil return kit	S225-05990-92




Fig. 8-4

■ Connecting the power supply cable

- 1 Connect the instrument's power cable to a power supply that is fitted with an earth leakage breaker.
- 2 Turn the instrument's power switch ON.




Fig. 8-5

■ Connection to the PC

The LCMS-8030/LCMS-8040 and computer are connected with a USB cable.

- 1 Connect the USB socket on the rear face of the LCMS-8030/LCMS-8040 and a USB socket on the computer.




Fig. 8-6

■ Connecting to the LC

When connecting to the CBM-20A or CBM-20Alite system controller:

- 1 Connect the start cable between the EVENT IN1 terminal on the rear face of the LCMS-8030/LCMS-8040 and the OUT1 or OUT2 terminal on the rear face of the CBM-20A/20Alite.




Fig. 8-7

8.1.3 About the Gas Used

This instrument uses a large volume of nitrogen gas. This poses a risk of anoxia symptoms (oxygen deficiency). Read this section carefully and give adequate consideration to the use of nitrogen gas.

■ Preventing anoxia symptoms (oxygen deficiency)

In accordance with the regulations for preventing anoxia symptoms, provide ventilation to maintain the concentration of oxygen at least 18 % at the installation site. (The normal concentration of oxygen in the air is about 21 %).

Install the instrument in a room that has a ventilation mechanism such as the type of draft chamber in general use (approx. 20 m³/min), and feed the exhaust tube into the draft chamber.

When using nitrogen gas, appoint a person in charge of operations where there is a danger of anoxia symptoms (a person who has completed a skills course for a first class superintendent of work involving a danger of anoxia symptoms) and give the workers special training relating to anoxia hazards.

To cover all eventualities, you are also recommended to use a separate oxygen analyzer. Details of the recommended portable oxygen analyzer are given below.

Recommended product: Type XO-2000 oxygen analyzer

- New Cosmos Electric Co., Ltd.
2-5-4 Mitsuyanaka, Yodogawa-ku, Osaka 532-0036, Japan
Tel: +81-6-6308-3111
FAX: +81-6-6308-8129
URL: <http://www.new-cosmos.co.jp/en/index.html>

■ Specifications of the nitrogen gas

In order to maintain the performance of the instrument, use gas with the specifications indicated below.

Nitrogen Supply pressure: 690 to 800 kPa

Purity: 97 % or greater

The maximum flow rate during use is 25 L/min.

You are recommended to use a nitrogen gas generator.

■ Specifications of argon gas

Argon gas is used for the CID gas. Use gas with the specifications indicated below.

Argon Supply pressure: 500 kPa or greater

Purity: 99.99 % or greater

**WARNING****Instructions**

- **Cautions on the use of high-pressure gas cylinders**
If you use a high-pressure gas cylinder as the gas supply source, follow the guidance of the cylinder dealership and others and handle the cylinders in a way that will avoid accidents.
 - High-pressure gas cylinders should be sited outdoors at a location where there is a good passage of air and no exposure to direct sunlight, and the gas should be fed indoors through tubing. In the case of liquefied gases in particular, these arrangements are obligatory under law.
 - Take care to ensure that high-pressure gas cylinders never reach high temperatures of over 40 °C. Also ensure that there will be no naked flames within 2 m of the cylinders.
 - Pay adequate attention to ventilation, and, as a start-of-work inspection, check for gas leakage using, for example, the soapy water test.
 - Secure high-pressure gas cylinders, e.g. with rope, so that they cannot be made to fall over/fall off.
 - Make sure that you use a pressure reducing valve of the oil prohibition type. In addition, do not use a valve that has oil adhering to the inside faces of the pipes that gas comes into contact with.
 - When the use of gas has finished, turn off the main cock on the high-pressure gas cylinder immediately.
 - Inspect the pressure gauge to check that it functions correctly at least once every three months.
 - Note that permission is required under the law to store more than 300 m³ (standard status) of high-pressure gases like these.
- Handle high-pressure gas cylinders correctly in compliance with regulations, the general high pressure gas preservation regulations and laws applicable in your country.

■ Supplying nebulizer gas and drying gas

- 1 Connect the nitrogen gas tubing to the rear of the instrument.
In the standard specifications, a $\phi 6$ mm tube joint is fitted. If, for example, you are connecting inch size tubing, prepare a convertor joint.
Use nitrogen gas with a purity of at least 97 %.




Fig. 8-8



NOTE

Connection to and disconnection from the one-touch joint

- To connect, insert the tube as far as it will go into the joint.
- To disconnect, press in the release bush.




Fig. 8-9

- 2 Open the main cock for the supply of nebulizer gas and drying gas.
Check that the supply pressure is 690 to 800 kPa.

■ Supplying CID gas

- 1 Connect the argon gas tubing to the rear of the instrument.
A Shimadzu M-shaped joint (see Fig. 8-8) is attached to the rear of the instrument.
The following gas conductor has been prepared for connecting the argon gas cylinder to the LCMS-8030/LCMS-8040.

Part Name	Part No.
Gas conductor (length: 2.5 m)	S201-48067

- 2 Open the argon gas main cock in order to supply CID gas.
Check that the supply pressure is 500 kPa.

About tubing joints

The parts described below are options. Purchase them if necessary.

■ Conversion from inch sizes

Part Name	Part No.
Tube joint 6.4	S225-15849-91




Fig. 8-10

■ Carrier gas tubing: Conversion from a Shimadzu M type joint

Part Name	Part No.
Gas pipe adapter kit	S225-10766-91




Fig. 8-11

8.1.4 Tubing for the LC Flow Line

■ Piping example

A typical example of the tubing used with the LCMS-8030/LCMS-8040 is shown below.




Fig. 8-12

■ For ESI (DUIS)




Fig. 8-13

■ For APCI




Fig. 8-14

8.1.5 Installing the Standard Sample

The standard sample is used for optimization of the supply pressure to the lens system and quadrupole rods inside the instrument in auto tuning, for sensitivity adjustment, for resolution adjustment, and for mass calibration.

CAUTION	
	CAUTION
	Turn the standard sample OFF from LabSolutions before removing tubing or the bottle cap.
Instructions	

- 1 Pour 40 to 80 mL of standard sample into the standard sample bottle.



Reference

["7.6.1 Replacing the Standard Sample" P.154](#)

- 2 Open the front door and set the standard sample bottle in the instrument.





Fig. 8-15



NOTE

A glass capillary tube is inserted into the resistance tube.

If bent abruptly, it could break.

Maintain a moderate bending radius of at least 40 mm while handling the tube.

- 3 Connect the gas tubing for delivering the standard sample.

4 Connect the resistance tube and the tubing to the ionization probe.




Fig. 8-16



NOTE

Connection to and disconnection from the one-touch joint

- To connect, insert the tube as far as it will go into the joint.
- To disconnect, press in the release bush.




Fig. 8-17



Reference

For a figure showing how to run the tube to the instrument body, see [Fig. 3-53](#).

8.1.6 Installing the Waste Container

! WARNING



For details on handling the waste container, be sure to read "[Static Electricity Precautions](#)" P.xiii.

Instructions

! CAUTION



- Install the waste container at a position lower than the surface on which the LCMS-8030/LCMS-8040 is installed.

Instructions

The waste liquid is expelled through the mechanism of a head drop. This will cause the waste liquid to overflow into the interior of the instrument.

- Make sure that the waste tube is not submerged below the surface of the liquid. This will cause the waste liquid to overflow into the interior of the instrument.

- 1 Prepare a waste container of an appropriate capacity (diameter of mouth: 29 to 31 mm) and place this container at a position lower than the instrument (e.g. on the floor).

Connect the drain in the manner shown [Fig. 8-18](#).




Fig. 8-18

8.2 Manual Tuning

In manual tuning, the instrument can be adjusted manually while monitoring peaks. If peaks cannot be checked while using manual tuning, execute auto tuning instead.



NOTE

- Ion monitoring using the standard sample can be accomplished in ESI analysis.
- In the APCI mode the flow volume is low and ions cannot be monitored stably.



CAUTION



Create a new tuning file when replacing the detector and never use a previous tuning file.

Instructions

8.2.1 Checking Ions in the [MS Tuning] Window

1

Connect the standard sample.



Reference

["8.1.5 Installing the Standard Sample" P.198](#)

2


Start up LabSolutions.



Reference

["3.1.2 Starting the Vacuum System" P.42](#)


- 3 Click the  (Tuning) icon.
The [Tuning] window is displayed.



The screenshot shows the 'Realtime Analysis' software interface. The 'Tuning' window is active, displaying three chromatograms: 'LCReady MSReady', 'MS Running Time: 0.00 / 10.00 min Scan#: 0 Inten.: 0', and 'Event# 1 Polarity: + Mode: MRM'. Below the chromatograms, the 'Instrument Parameters View' is shown, including 'MS Simple Settings' and 'AutoPurge' options. The 'End Time' is set to 10.000 min. The 'Event# 1' table shows a single entry for MRM at 0.000-10.000 min.

Item	Value	Setting	Units
Interface	ESI		
Nebulizing Gas Flow	1.5	1.5	L/min
Drying Gas Flow	15.0	15.0	L/min
Detector Voltage		1.50	kV
DL Temperature	250	250	C
Heat Block Temperature	400	400	C
Interface Voltage		4.5	kV
Interface Current	0.1		uA
PG Vacuum	1.1e+002		Pa
IG Vacuum	3.0e+004		Pa
CID Gas	17	17	kPa
Mode	Isocratic flo	Isocratic flow	
Total Flow	0.200		mL/min
Flow Rate	0.000		L
Pump A Flow	0.050	0.000	mL/min
Pump B Flow	0.150	0.000	mL/min
Injection Volume			uL



- 4 Click the  (Manual Tuning) icon.
The [MS Tuning] window is displayed.



The screenshot shows the 'Realtime Analysis' software interface with the 'MS Tuning' window active. The window displays four intensity vs. m/z plots for different ions. The 'Tuning Condition' is set to 'Manual Tuning'. The 'Tuning Mode' is also set to 'Manual Tuning'.

m/z	Inten.	FWHM									
168.10	0	0.00	256.15	0	0.00	344.20	0	0.00	652.40	0	0.00
Mag.: 20			Mag.: 20			Mag.: 20			Mag.: 20		

- 5 Turn ON the nebulizer gas, the heaters (DL, heated block) and the [MS detector] button.



- 6 Turn ON the  (standard sample) button. The sample will not flow out for approximately one minute while the sample tubing becomes filled with the sample.

7 Check the peaks in the [MS Tuning] window.

In manual tuning, screen settings are configured so that the following ions can be monitored in each analysis mode.

- Ions configured in positive mode

Analysis Mode	Peak Monitor Settings m/z
Q1 and Q3 scan	168.10, 256.15, 344.20, 652.40, 1004.60, 1224.75, 1603.15, 1893.40
Precursor ion scan	168.10 > 89.00, 300.20 > 89.00, 344.20 > 177.10, 388.25 > 177.10
Product ion scan	168.10 > 45.05, 168.10 > 89.00, 388.25 > 133.10, 388.25 > 177.10
Neutral loss scan	256.15, 300.20, 344.20, 388.25, neutral loss 211.15

- Ions configured in negative mode


Analysis Mode	Peak Monitor Settings m/z
Q1 and Q3 scan	503.15, 1007.30
Precursor ion scan	503.15 > 89.00, 1007.30 > 503.15
Product ion scan	503.15 > 89.00, 503.15 > 179.05
Neutral loss scan	503.15 neutral loss 414.15, 503.15 neutral loss 324.10



Reference


["8.3.3 Mass of the Standard Sample" P.215](#)

Switching the number of windows displayed



Item	Value	Setting	Units
Interface	ESI		
Nebulizing Gas Flow	1.5	1.5	L/min
Drying Gas Flow	15.0	15.0	L/min
Detector Voltage	2.32		kV
DL Temperature	250	250	C
Heat Block Temperature	400	400	C
Interface Voltage	4.5		kV
Interface Current	3.2		uA
PG Vacuum	9.7e+001		Pa
IG Vacuum	1.1e-003		Pa
CID Gas	230	230	kPa
Mode	isocratic flo	isocratic flow	
Total Flow	0.200		mL/min
B Conc.	0.0		L/7500
Pump A Flow	0.050	0.000	mL/min
Pump B Flow	0.150	0.000	mL/min
Injection Volume			uL

Switching between the positive ion and negative ion monitor



Item	Value	Setting	Units
Interface	ESI		
Nebulizing Gas Flow	1.5	1.5	L/min
Drying Gas Flow	---	15.0	L/min
Detector Voltage	2.32		kV
DL Temperature	250	250	C
Heat Block Temperature	400	400	C
Interface Voltage	4.5		kV
Interface Current	0.9		uA
PG Vacuum	9.7e+001		Pa
IG Vacuum	1.0e-003		Pa
CID Gas	230	230	kPa
Mode	isocratic flo	isocratic flow	
Total Flow	0.200		mL/min
B Conc.	0.0		L/7500
Pump A Flow	0.050	0.000	mL/min
Pump B Flow	0.150	0.000	mL/min
Injection Volume			uL

■ Entering the m/z value to be monitored

The screenshot shows the 'Realtime Analysis' software interface. The main window displays four mass spectra plots. The first plot is selected, and its m/z value (168.10) is entered into the 'm/z:' field. The 'Mag.' field is set to 'ZZ'. The software also displays various instrument parameters and a 'Details...' table on the right.

Item	Value	Setting	Units
Interface	ESI		
Nebulizing Gas Flow	1.5	1.5	L/min
Drying Gas Flow	---	15.0	L/min
Detector Voltage	---	2.32	kV
DL Temperature	250	250	C
Heat Block Temperature	400	400	C
Interface Voltage	---	4.5	kV
Interface Current	0.9		uA
PG Vacuum	9.7e+001		Pa
IG Vacuum	1.0e+003		Pa
CID Gas	230	230	kPa
Mode	Isocratic flo	Isocratic flow	
Total Flow	0.200		mL/min
B-Comp	0.0		
Pump A Flow	0.050	0.000	mL/min
Pump B Flow	0.150	0.000	mL/min
Injection Volume			uL

■ To change the monitoring magnification



Enter a value at [Mag.].

The screenshot shows the 'Realtime Analysis' software interface. The main window displays four mass spectra plots. The first plot is selected, and its 'Mag.' field is now set to '22', indicating a change in magnification. The rest of the interface remains the same.

Move the mouse cursor into the ion observation panel and double click to display the data standardized at the maximum intensity.


8.2.2 Exiting the [MS Tuning] Window

- 1 Click the  (Result Tuning) icon to exit manual tuning.



The save confirmation message is displayed. Click [Yes] to save the result of manual tuning.

- 2 Select whether you want to save the peak profile and spectrum of the standard sample observed in manual tuning or not. If [Yes] is selected, the following message will be displayed.




- 3 Select whether the saved tuning file is to be set as the default tuning file or not.



Reference

"3.8.2 Starting Auto Tuning" P.82



8.2.3 Explanations of Parameters




Fig. 8-19

■ Gas

Parameter Name	Default Value	Maximum	Minimum	Remarks
Nebulizer gas	3.0 L/min (ESI) 2.0 L/min (DUIS) 3.0 L/min (APCI)	3.0 L/min (ESI, DUIS) 4.4 L/min (APCI)	0.5 L/min (ESI, DUIS) 0.5 L/min (APCI)	
Drying gas	15 L/min (ESI, DUIS) 5 L/min (APCI)	20 L/min	3 L/min	Execute with auto tuning OFF
CID gas	230	230	17	0 when turned OFF in the [System Control] window

■ High voltage

Parameter Name	Default Value	Maximum	Minimum	Remarks
Interface HV	ESI (+) 4.5 kV (-) -3.5 kV APCI (+) 4.5 kV (-) -3.5 kV	(+) 5.0 kV (-) -5.0 kV	(+) 0.0 kV (-) 0.0 kV	
DUIS corona HV	(+) 4.5 kV (-) -3.5 kV	(+) 5.0 kV (-) 5.0 kV	(+) 0.0 kV (-) 0.0 kV	

■ Heater

Parameter Name	Default Value	Maximum	Minimum	Remarks
DL heater	250 °C	300 °C	50 °C	
Heated block	200 °C (APCI) 400 °C (ESI, DUIS)	300 °C (APCI) 500 °C (ESI, DUIS)	50 °C	
APCI heater	350 °C	500 °C	50 °C	

■ Lens voltage

Parameter Name	Default Value	Maximum	Minimum	Remarks
DL	0 V	(+) 100 V (-) -100 V	0 V	
Qarray DC	0 V	100 V	-100 V	
Qarray RF	1) Depends on m/z	130 V	10 V	
Skimmer	0 V	0 V	0 V	Fixed
Multipole 1 DC	0 V	0 V	0 V	
Multipole 2 DC	(+) -1 V (-) +1 V	10 V	-10 V	
Multipole RF	2) Depends on m/z	300 V	10 V	Same for multipole 1 and 2
Entrance lens 1	(+) -1 V (-) +1 V	10 V	-10 V	
Entrance lens 2	(+) -2 V (-) +2 V	50 V	-50 V	

1) Qarray RF




Fig. 8-20

2) Multipole RF




Fig. 8-21

■ Analysis unit

Parameter Name	Default Value	Maximum	Minimum	Remarks
Q1 pre-rod DC	1) Depends on m/z	50 V	- 50 V	
Q1 main rod DC	(+) - 5 V (-) + 5 V	50 V	- 50 V	
Q1 post-rod DC	(+) - 5 V (-) + 5 V	50 V	- 50 V	
Collision cell RF (CCRF)	2) Depends on m/z	300 V	0 V	Optimized for each analysis mode
Collision energy	(+) - 15 V (-) 30 V	180 V	- 180 V	
Q3 pre-rod DC	3) Depends on m/z	50 V	- 50 V	Difference with CE
Q3 main rod DC	4) Depends on m/z	50 V	- 50 V	Difference with CE

1) Q1 pre-rod DC with positive ions

With negative ions, the positive and negative voltages are interchanged.




Fig. 8-22

2) Collision cell RF




Fig. 8-23

3) Q3 pre-rod DC with positive ions

With negative ions, the positive and negative voltages are interchanged.




Fig. 8-24

4) Q3 main rod DC with positive ions

With negative ions, the positive and negative voltages are interchanged.




Fig. 8-25

■ Detector

Parameter Name	Default Value	Maximum	Minimum	Remarks
Conversion dynode	(+) -6 kV (-) 6 kV	(+) -10.0 kV (-) 10.0 kV	0 V	
Detector	-1.75 kV	-3.5 kV	0 V	

8.3 Standard Sample

For carrying out sensitivity adjustment and atomic mass number adjustment of the instrument, the LCMS-8030/LCMS-8040 uses a mixture of polyethylene glycol (PEG) 200, 600, 1000, 2000 polypropylene glycol (PPG) and raffinose solution as the standard sample.

Sample for instrument calibration P/N S225-14122-01 (200 mL)

Using this standard sample, auto tuning automatically performs sensitivity optimization and mass calibration.

8.3.1 Method for Preparing the Standard Sample

The method for preparing the standard sample is as follows.

The result is the same as P/N S225-14122-01 (200 mL).

1 Make up the dilution solvent.

Dilution solvent (when making approx. 1 L)

Pure water	800 mL
Methanol	200 mL
Ammonium acetate	14.5 mg

2 Make up the stock solution.

Dissolve the sample compounds indicated below in 100 mL of the dilution solvent described in 1 above.

PEG200	0.75 µL
PEG600	1.0 µL
PEG1000	150 µL
PPG2000	100 µL

This is the stock solution.

(It is 100 times the concentration of the standard sample. However, it does not contain raffinose.)

PEG1000 is a solid at room temperature. To weigh it, warm it to around 60 °C to liquefy it, weigh it quickly in a disposable micro pipette and then dissolve it in the solvent.

3 Make up the standard sample for auto tuning.

- 1 Dilute the stock solution created in step 2 to 1/100th concentration using the dilution solution created in step 1.
- 2 Dissolve raffinose in this solution to achieve a concentration of 15 mg/L.

This completes the procedure.

The resulting concentrations are as indicated below.

PEG 200	: 0.075 µL/L
PEG 600	: 0.1 µL/L
PEG 1000	: 15.0 µL/L
PPG 2000	: 10.0 µL/L
Raffinose	: 15 mg/L

NOTE

The PEG and PPG used for the standard sample easily stick to the equipment and are difficult to remove even by washing, so the pipette, containers, etc. used for weighing these compounds should not be used for the preparation of the mobile phase or reagents.

This will cause contamination of reagents and background noise during analysis.

8.3.2 Reagents

PEG200	500 g	Wako Pure Chemical Industries	P/N 167-09045
PEG600	500 g	Wako Pure Chemical Industries	P/N 168-09075
PEG1000	500 g	Wako Pure Chemical Industries	P/N 165-09085
PPG2000	500 g	Wako Pure Chemical Industries	P/N 164-05895
Raffinose	5 g	Wako Pure Chemical Industries	P/N 182-00011
Solvent			
Methanol for LC use	1 L		P/N 017-40019-06
Distilled water for LC use	1 L		P/N 017-40513-01
Ammonium acetate, special grade	500 g	Wako Pure Chemical Industries	P/N 019-02835

8.3.3 Mass of the Standard Sample

For the positive ionization mode

PEG (Polyethylene Glycol): $\text{HOCH}_2(\text{CH}_2\text{OCH}_2)_n\text{CH}_2\text{OH}$

Auto tuning is performed by using the PEG + NH_4 molecules with added ammonium ions.

n=1	2	3	4	5	6	7	8	9
124.0974	168.1236	212.1498	256.176	300.2022	344.2284	388.2547	432.2809	476.3071
10	11	12	13	14	15	16	17	18
520.3333	564.3595	608.3857	652.4119	696.4382	740.4644	784.4906	828.5168	872.543
19	20	21	22	23	24	25	26	27
916.5692	960.5954	1004.622	1048.648	1092.674	1136.7	1180.727	1224.753	1268.779




Fig. 8-26

For the negative ionization mode

Raffinose: $\text{C}_{18}\text{H}_{32}\text{O}_{16}$

Auto tuning performs mass calibration using Raffinose 503.16 (M-H) deprotonated molecules and 1007.33 (2M-H) dimer ions.





Fig. 8-27

8.4 Example Configuration of a High-Pressure Gradient System (Auto Injector)

■ Prominence system

Flow of solvent



Function of each unit




Fig. 8-28

9

Specifications

9.1 Performance

Item	Specifications
Mass range	m/z 10 to 2000
Resolution	R = 2 M (however, 50 % valley, 1000 u/sec.)
Sensitivity (MRM)	ESI Positive ions Reserpine 1 pg, m/z 609.3 > 195 S/N > 200 (RMS) @LCMS-8030 S/N > 1000 (RMS) @LCMS-8040
Maximum scan speed	15000 u/sec
Positive/negative switching time	15 msec
MRM settings	Minimum Dwell Time: 1 msec Minimum Pause Time: 1 msec

9.2 Hardware

Item	Specifications
Interface section	ESI probe / APCI probe (option) / DUIS (option) and ion sampling unit Drying gas Max. flow rate 20 L/min
ESI probe	Sample flow rate (max.) 2 mL/min Probe voltage ± 5 kV max. Nebulizer gas Flow rate 3.0 L/min max.
APCI probe (option)	Sample flow rate (max.) 2 mL/min Corona needle voltage ± 5 kV max. Probe temperature 500 °C max. Nebulizer gas Flow rate 4.4 L/min max.
DUIS (option)	Sample flow rate (max.) 1 mL/min Probe voltage ± 5 kV max. Corona needle voltage ± 5 kV max.
DL	Temperature 300 °C max. Voltage ± 100 V
Heated block	Temperature 500 °C max. (ESI, DUIS) 300 °C max. (APCI)
Analysis/detection section	
Analysis rods	Hyperbolic quadrupole rods made of molybdenum, with pre-rods
Collision cell	Multipole type ultra high speed collision cell
Detector	Secondary electron multiplier with conversion dynode
Vacuum system	Triple inlet turbo molecular pump, 1 unit (triple inlet) Vacuum speed 40 L/s + 260 L/s + 210 L/s Rotary pump, 1 unit Vacuum speed 28 m ³ /h

9.3 Software

9.3.1 Analysis

Item	Specifications
Analysis conditions	Batch register the measurement conditions for the LC, the measurement conditions for the MS, and the conditions for data processing in the Method, and register the report output format in the template file.
MS analysis	Q3 scan (max. 512 events), Q3 SIM (max. 512 events × 32 channels) Q1 scan (max. 512 events), Q1 SIM (max. 512 event × 32 channels)
MS/MS analysis	MRM (max. 512 events × 32 channels) Product ion scan (max. 512 events) Precursor ion scan (max. 512 events) Neutral loss scan (max. 512 events)
Automatic analysis	Automatic analysis can be performed through the combination of the auto injector (option) and the schedule function.

9.3.2 Data Processing

Item	Specifications
Data processing	Drawing of LC chromatograms, MS chromatograms and MS spectra Elimination of background noise in LC chromatograms Background subtraction and equalization for MS spectra Area calculations for LC chromatograms and MS chromatograms Calculation of column performance from LC chromatograms
Quantitative calculations	For both LC data and MS data Number of identified peaks: Max. 3000 Identification method: Absolute retention time method / relative retention time method, time band method /time window method, reference ions can be used (maximum of 5). Quantitative calculations: Correction percentage method (possible to add a scale factor), internal standard method, absolute calibration curve method (external standard), standard addition method Calibration curve: Straight line (method of least square, averaging method), broken line, secondary curve, tertiary curve, maximum of 64 calibration points (averaging up to 10 times is possible for each calibration point), weighting
Library search	Database: Private library Search modes: Similarity search, index search Number of libraries as simultaneous search targets: Max. 5
Library editing	Editing of the private library
Report output	LC chromatograms, MS chromatograms, MS spectra, profiles, LC peak tables, MS peak tables, MS status log, LC column performance reports, library search results, summary reports
Batch processing	Data processing by the scheduler, and continuous automatic processing of report output are possible.

9.3.3 Instrument Control

Item	Specifications
Instrument Control	Starting and stopping the MS unit Automatic and manual adjustment of the MS unit LC unit and MS unit diagnosis functions

9.3.4 User Management

Item	Specifications
User Management	A maximum of 37 levels of privileges can be assigned to each user.

9.4 Installation Conditions

Item	Specifications
Required power supply	<p>Single phase 230 V AC (50 Hz/60 Hz), 15 A, for LCMS-8030/LCMS-8040</p> <p>Power for LC</p> <p>Refer to the specifications or the instruction manual provided with your LC. (The current capacity varies according to the configuration of the LC unit.)</p> <p>*) For details on the power supply specifications of the computer, printer and peripheral equipment, check either the specifications or the instruction manuals of the devices you are using. The specifications given below are a rough guide.</p> <p>Power for PC</p> <p>Supply correct voltage to the PC or printer, etc. referring to the instruction manuals accompanying the PC or printer, etc.</p> <div style="border: 1px solid black; padding: 5px; margin-top: 10px;"> <p style="text-align: center;"> CAUTION</p> <p style="text-align: center;"> Use a PC, printer, etc. that conform to your local laws and regulations.</p> <p>Instructions</p> </div> <p>*) As far as voltage fluctuation is concerned, the instrument is designed to operate adequately well provided the range of voltage fluctuation is within $\pm 10\%$ including the amplitude of the fast noise superimposed in the AC line, but the voltage fluctuation range required in order to guarantee the performance specifications is $\pm 5\%$ including noise.</p> <p>Grounding: 100 Ω max.</p>
Installation environment	<p>Temperature 18 °C to 28 °C</p> <p>Humidity 40 % to 70 % (to be no dew condensation)</p> <p>Other Apart from this, the instrument should be installed in an environment with few disturbances such as dust, vibration, electromagnetic wave noise, corrosive gases, interfering magnetic fields and so on.</p>
Exhaust gas	<p>In order to process the exhaust gas of the rotary pump safely, an exhaust system such as a draft chamber is required.</p>
Required gas	<p>Nitrogen gas</p> <p>Supply pressure: 690 to 800 kPa</p> <p>Purity: 97 % or greater</p> <p>Max. flow rate during use: 25 L/min</p> <p>Argon gas</p> <p>Supply pressure: 500 to 800 kPa</p> <p>Purity: 99.99 % or greater</p> <p>Max. flow rate during use: Approx. 10 mL/min</p>

9.5 Example Installation

	Dimensions			Weight, kg
	W	D	H	
LCMS-8030 LCMS-8040	1178	530	561	130
Rotary pump	597	180	276	44




Fig. 9-1

10

Maintenance Parts

10.1 Consumables

Part Name	Part No.	Remarks
ESI probe ASSY		
CAPILLARY ASSY	S225-14948-91	Set of tubes for ESI
PEEK TUBE 5999	S228-32999-01	
ELBOW, KQ2L23-M5-ROHS	S035-60724-21	
APCI probe ASSY		
APCI PIPE ASSY	S225-15845-91	SUS capillary tube
HEATER ASSY	S225-15619-41	APCI heater unit (with ferrule)
FERRULE	S225-03748-03	
HALF UNION, KQ2H01-M5-ROHS	S035-60725-01	
Corona needle ASSY		
NEEDLE ASSY	S225-15877-92	Needle (for both APCI and DUIS)
Interface ASSY		
DL PIPE ASSY	S225-15718-91	DL unit set
ORIFICE	S225-15479	
SAMPLING CONE	S225-15487	
WASHER PEEK	S023-65106-01	Seal for fixed part of heater flange
SI unit ASSY		
TUBE, SI	S225-15848-91	Resistance tube, joints, tubing for LC
SPACER, FKM	S225-15697-01	Glass bottle packing
PEEK FRIT	S228-48607-91	Made of PEEK
Standard sample bottle	S038-00512-01	
Bottle cap	S225-15868-91	Standard sample bottle cap
Standard sample	S225-14122-01	
Detection unit		
EM, MS644	S225-14168-01	Secondary electron multiplier
Vacuum system		
RP Oil Ultragrade19 (4L)	S017-30163-02	About 1.5 L used at each oil change
IG GAUGE	S225-09490-01	Vacuum gauge
FILAMENT FOR PB1, ROHS	S225-20310-91	Vacuum gauge (Pirani gauge)

10.2 Spare Parts

Part Name	Part No.	Remarks
Waste tube		
SILICON RUBBER TUBE, 7X10NL	S016-31350-19	Piping from the leak tray
PVC TUBE, R3603 1/2X3/4X1/8	S016-31414	Drain tube from the spray unit, for use outside the instrument
TUBE, SE-200 1/2-3/4	S016-37619-02	Drain tube from the spray unit, for use inside the instrument
Vacuum system		
TMP, split flow 310	S225-14179-01	Triple inlet turbo molecular pump body + power supply
ROTARY PUMP, E2M28	S225-09309-02	For 230 V
CABLE ROTARY PUMP	S225-17224-41	Power supply cable for the rotary pump
LEAK VALVE ASSY	S225-03540-94	
HOSE TG-32	S018-31555-06	Use a 1.5 m RP intake hose (purchase unit: m)
HOSE, CHEMIFLEX 15MM	S016-31697-02	Use a 0.5 m TMP exhaust hose (purchase unit: m)
SLEEVE, PTFE 15X22 CL	S018-31511	Use a 5 m RP exhaust hose (purchase unit: m)
HOSE BAND, GEAR S 50	S037-61064	For spring hose #32
HOSE BAND, WIRE 24	S037-61023	For spring hose #15
HOSE BAND, WIRE SY-22	S037-61002	For vinyl hose 15 × 1.5
CENTER RING, KF10	S035-06004-21	
CENTER RING, KF25	S035-06004-24	
CLAMP, KF16	S035-06004-01	
CLAMP, KF25C	S035-06004-02	
FLANGE PIPE, TQ-VAC	S225-12211-91	For RP intake
ELBOW, KF25-#15	S225-03535-91	For TMP back (TMP side)
FLANGE	S202-55374	For RP intake side (RP side)
O-RING, 4D G30	S036-12502	For RP intake
Housing ASSY		
O-RING, 4D P15	S036-11216	
O-RING, 4D P40	S036-11243	
O-RING, 4D-S56	S036-19004-41	
O-RING, AS568A-253 4D	S036-15552-53	
O-RING, AS568A-272 4D	S036-15552-72	
O-RING, AS568A-278 4D	S036-15552-78	
O-RING, AS568A-341 4D	S036-15553-41	
O-RING, 4D P18	S036-11218	
O-RING, 4D P5	S036-11203	
PACKING 4TX14X22	S261-00207-02	For the ion gauge

Part Name	Part No.	Remarks
Interface section		
IF FLANGE	S225-15477-41	Interface flange body
INSULATOR FLANGE	S225-15481	
PLATE, ORIFICE	S225-15710	
SOCKET ASSY	S225-15717-41	
BOLT, SST HEXSOCKET M6X16	S022-27104	
O-RING, 4D P4	S036-11202	
O-RING, 4D S8	S036-19004-05	
O-RING, 4D S14	S036-19004-11	
O-RING, 4D-S70	S036-19004-46	
HEATER FLANGE ASSY	S225-15486-41	Heater flange body
BLOCK, HEATER	S225-15488	
INSULATOR BLOCK	S225-15489	
HEATER ASSY	S225-15716-41	
INSULATOR, FLANGE	S225-15491	
BUSH, HEATER	S225-15498	
O-RING, 4D S100	S036-19004-53	
O-RING, 4D S18	S036-19004-14	
O-RING, 4D S14	S036-19004-11	
O-RING, 4D S10	S036-19004-07	
Ion source unit		
DOOR ASSY	S225-12100-41	
WINDOW	S225-12106	Window glass
DOOR COVER	S225-12101	Plastic cover
O-RING, 4D G145	S036-19004-53	
O-RING, 4D P14	S036-11215	
O-RING, 4D G50	S036-12506	
O-RING, 4D G100	S036-12517	
O-RING, 4D P26	S036-11228	
ESI probe		
ESI PROBE ASSY	S225-14949-41	ESI probe, complete
CAPILLARY	S225-14915	
NOZZLE	S225-14902	
SUS COUPLING	S228-16447-03	
COUPLING	S225-14903	
MALE NUT, PEEK	S228-18565-84	5 included
HIGH VOLTAGE CABLE	S225-14947-41	
O-RING, 4D P7	S036-11205	
O-RING 4D P9	S036-11207	

Part Name	Part No.	Remarks
ESI probe		
O-RING, 4D-S8	S036-19004-05	
ELBOW, KQ2L23-M5-ROHS	S035-60724-21	
APCI probe		
APCI PROBE ASSY	S225-14271-41	APCI probe, complete (excluding corona needle)
JOINT ASSY	S225-15788-91	
NUT	S225-15739	
HIGH VOLTAGE CABLE	S225-15888-41	Common to APCI and DUIS
ADAPTER	S225-04993	
NEEDLE UNIT ASSY	S225-14290-41	Common to APCI and DUIS
HALF UNION, KQ2H01-M5-ROHS	S035-60725-01	Gas joint
DUIS		
NEEDLE UNIT ASSY	S225-14290-41	Common to APCI and DUIS
HIGH VOLTAGE CABLE	S225-15888-41	Common to APCI and DUIS
Analysis unit		
ROD ASSY	S225-03550-91	Quadrupole rod (with pre-rod)
PRE-ROD ASSY	S225-03560-91	Common to Q1 and Q2
PRE-ROD ASSY	S225-03560-92	Common to Q1 and Q2
SPRING ASSY	S225-12257-91	Common to Q1 and Q2
PRE-ROD	S225-03552-01	Q1 post-rod
INSULATOR	S225-03553	For post-rod (pre-rod) coupling
NUT	S225-03554	For post-rod (pre-rod) coupling
PIN, POST ROD	S225-12258-42	Post-rod short
Q1 CABLE ASSY	S225-12270-91	Q1 cable
Q3 CABLE ASSY	S225-12270-92	Q3 cable
BAND	S225-03566-03	Quadrupole rod tension plate
CID cell		
COLLISION CELL ASSY	S225-12070-41	CID cell (for the LCMS-8030)
	S225-12070-42	CID cell (for the LCMS-8040)
BAND CID CELL	S225-12015	CID cell tension plate
COLLISION CELL CABLE ASSY	S225-14150-91	
O-RING, 4D-S4	S036-19004-01	
COLLISION CELL Lens1	S225-12217-41	
COLLISION CELL Lens2	S225-12218-41	
COLLISION CELL Lens3	S225-12219-41	
COLLISION CELL Lens out	S225-12220-41	
COLLISION CELL Plate	S225-14092	

Part Name	Part No.	Remarks
Detection unit		
Detector ASSY	S225-14168	Detector, complete
EM, MS644	S225-14168-01	Secondary electron multiplier only
EM HV CABLE ASSY	S225-12292-91	Cable for detector voltage supply
EM SIG CABLE ASSY	S225-12293-91	Signal cable
Lens system		
Qarray ASSY	S225-12030-91	Qarray + skimmer
SKIMMER	S225-12034-01	
OP1 ASSY	S225-12040-91	Qarray side multipole ASSY (for the LCMS-8030)
QP1 ASSY	S225-13646-41	Qarray side multipole ASSY (for the LCMS-8040)
OP2 ASSY	S225-12050-91	Q1 side multipole ASSY (for the LCMS-8030)
QP2 ASSY	S225-13648-41	Q1 side multipole ASSY (for the LCMS-8040)
OP1 Lens	S225-12045-01	Entrance lens of the QA side multipole ASSY (for the LCMS-8030)
QP1 Lens	S225-13630-01	Entrance lens of the QA side multipole ASSY (for the LCMS-8040)
OP2/QP2 Lens	S225-12052-01	Entrance lens of the Q1 side multipole ASSY (for the LCMS-8030/8040)
SI unit		
MANIFOLD, SI	S225-15695-91	Manifold
HALF UNION, GWJS6-M5	S035-65415-06	Gas joint
HOUSING, LINEFILTER	S228-46358	
PEEK TUBE 1.6X0.25	S228-32999-03	Use 95 mm (purchase unit: m)
MALE NUT, PEEK	S228-18565-84	5 included
BRACKET, SI	S225-15560	Mounting plate
Gas control system		
GAS Controller ASSY	S225-12008-41	Gas control unit set
FEP TUBE 1/16, NEB	S225-14255-41	For nebulizer gas
FEP TUBE 4 mm, DRY	S225-14255-42	For drying gas
FEP TUBE 1/16, DRY	S225-14255-43	For drying gas
POLYURETHANE TUBE, U2-4-6X4BK	S016-46021	
VALVE ASSY (NEB)	S225-51226-01	
VALVE ASSY (DRY)	S225-14292-41	
VALVE ASSY SAGINOMIYA	S221-48813-91	
FEP TUBE 6mm, GAS	S225-15846-93	
TUBE JOINT, M	S225-10766-91	For conversion from M type (option)
TUBE JOIN, 6.4	S225-15849-91	For conversion from 6.4 mm (option)
FILTER, 2300B-SS-1/8-2U	-	
REGULATOR, AR10-M5BG	S040-72549-51	For CID gas pressure reduction
PCB ASSY, TQ-FLOW	S225-17010-42	PCB (for maintenance, adjusted)

Part Name	Part No.	Remarks
CID gas vacuum introduction unit		
CID VALVE ASSY	S225-12103-41	
CID GAS INLET ASSY	S225-12131-41	CID gas introduction unit, complete
CID GAS CAPILLARY	S225-12131-42	Capillary and connection parts only
Fan		
FAN-ASSY	S225-14263-41	Main power supply, TMP cooling fan
PROBE FAN ASSY	S225-14190-41	Fan for inside the probe cover
FAN60MM ASSY	S225-09853-41	Fan for CID-RF power supply
RF FAN ASSY	S225-12151-41	For the main RF power supply
FILTER	S042-60935-14	Filter media only
Cover		
TOP COVER	-	Top panel
FRONT COVER	-	Cover for entire surface
SIDE PANEL R	-	Right side panel
SIDE PANEL L	-	Left side panel
PROBE COVER ASSY	-	Probe cover, complete
TQ PROBE COVER FRONT	S225-12233	Plastic cover for the probe cover front section only
GAS SPRING	S034-33018-05	Gas spring for the probe cover section
OP COVER	-	Multipole maintenance section cover
SI_DOOR	S225-12236	SI cover (for the LCMS-8030)
	S225-12236-01	SI cover (for the LCMS-8040)

■ Electrical parts

Part Name	Part No.	Remarks
Main power supply		
PCB ASSY, TQ-HPSCONT	S225-16940-40	
PCB ASSY, TQ-DCDCPS	S225-16950-40	
POWER SUPPLY, ZWS240PAF-24/J H	S074-80435-38	
POWER SUPPLY, ZWS150AF-24/J HFP	S074-80429-75	
POWER SUPPLY, ZWS100AF-24/J	S074-80429-55	
POWER SUPPLY, ZWS30-5/J HFP	S074-80383-22	
HEATER TRANS ASSY	S225-13518-41	Heater transformer
FILTER, MAS-1215-33	S075-00073-13	Noise filter
PROTECTOR, BAM2151310	S065-90768-12	Power switch
CABLE, MAIN POWER	S225-17112-41	
Control system		
PCB ASSY, TQ-CPU	S225-16900-41	CPU
PCB ASSY, TQ-LENSDC	S225-16930-41	LENS
PCB ASSY, TQ-ANALOG	S225-16920-41	ANALOG
CDL transformer		
LCMS CDLTRANS2	S225-13515	For the CDL heater
High-voltage power supply		
TQ-HV SIG ASSY	S225-17099-41	High-voltage power supply unit for detector
HV IF SUB ASSY	S225-13347-41	High-voltage power supply unit for interface
CABLE, PR HV P OUTPUT	S225-13602-42	HV cable for IF
HV DUAL SUB ASSY	S225-13348-41	High-voltage power supply unit for DUIS
CABLE, NED HV N OUTPUT	S225-13603-42	HV cable for DUIS
CABLE, NED HV INPUT	S225-13599-41	Flat cable for DUIS
RF power supply		
TQ-MAIN-RF	S225-12281-42	Main RF power supply shared by Q1 and Q3
TQ-QP-OP PS ASSY	S225-12161-41	Qarray, multipole RF power supply
TQ-CID PS ASSY	S225-12162-41	CID cell RF power supply
PCB ASSY, TQ-CNOR	S225-17060-41	Relay PCB for Qarray, multipole RF power supply (for the LCMS-8030)
PCB ASSY, TQ-CNOR2	S225-17060-42	Relay PCB for Qarray, multipole RF power supply (for the LCMS-8040)
PCB ASSY, TQ-CNCR	S225-17070-41	Relay PCB for collision cell RF power supply (for the LCMS-8030)
PCB ASSY, TQ-CNCR2	S225-17070-42	Relay PCB for collision cell RF power supply (for the LCMS-8040)

Part Name	Part No.	Remarks
Detector		
BOARD, PRE AMP	S225-16960-41	
LED		
PCB ASSY, LCMS-LED	S225-13480-41	PCB ASSY, LCMS-LED
Start signal		
CABLE, LCMS EVENT	S225-17126-41	MS-LC event cable
Analog output		
CABLE, TQ-ANALOG	S225-17128-41	MS-FCV signal cable
USB cable		
CABLE, PC-MS ASSY	S225-17125-41	USB cable
Cable for FCV		
CABLE, TQ-FCV	S225-17127-41	MS-FRC analog cable

10.3 Unit Configurations

Part Name	Part No.	Remarks
ESI PROBE ASSY	S225-14949-41	
APCI PROBE	S225-14271-41	
NEEDLE UNIT ASSY	S225-14290-41	Common to APCI and DUIS
IF FLANGE ASSY	-	
HEATER FLANGE ASSY	S225-15486-41	
IF FLANGE	S225-15477-41	
DL PIPE ASSY	S225-15718-91	
ORIFICE PART	-	
SPRAY UNIT	-	
LENS UNIT	-	
DETECTOR UNIT	-	
Qarray ASSY	S225-12030-91	
OP1 ASSY	S225-12040-91	(for the LCMS-8030)
QP1 ASSY	S225-13646-41	(for the LCMS-8040)
OP2 ASSY	S225-12050-91	(for the LCMS-8030)
QP2 ASSY	S225-13648-41	(for the LCMS-8040)
VACUUM UNIT	-	
ION GAGE	-	
PIRANI GAGE	-	
DRAIN TUBE	-	
STANDARD SAMPLE SOURCE	-	
NEBULIZER GAS	-	
ACCESORY	-	
STARTUP KIT	S225-13915-42	
OIL RETURN KIT	S225-05990-92	

10.3.1 ESI Probe: S225-14949-41

No.	Part Name	Part No.	Remarks
1	CAPILLARY ASSY	S225-14948-91	Capillary ASSY consumables
1-1	CAPILLARY	S225-14915	Capillary tube body
2	NOZZLE	S225-14902	
3	COUPLING, 1.6	S228-16447-03	SUS
4	COUPLING	S225-14903	PEEK (for securing capillary tube)
5	HIGH VOLTAGE CABLE	S225-14947-41	
6	ELBOW, KQ2L23-M5-ROHS	S035-60724-21	
7	PEEK TUBE 5999	S228-32999-01	Use 420 mm (purchase unit: m)
8	MALE NUT, PEEK	S228-18565-84	5 included
9	O-RING, 4D P7	S036-11205	
10	O-RING, 4D S8	S036-19004-05	
11	SCREW, SST FLAT HEAD M3X16	S020-12109	For securing the nozzle
12	SET SCREW, SST M3X3 CUP POINT	-	For securing the high-voltage cable
13	NUT, KNURL	S225-14914	For securing the capillary tube
14	NUT, RID-6A	S228-16034	For securing the SUS, coupling




Fig. 10-1

10.3.2 APCI Probe: S225-14271-41

No.	Part Name	Part No.	Remarks
1	APCI PIPE ASSY	S225-15845-91	
2	JOINT ASSY	S225-15788-91	
3	HEATER ASSY	S225-15619-41	
4	ADAPTER	S225-04993	
5	FERRULE	S225-03748-03	
6	NUT	S225-15739	
7	SCREW, SST SEMS P3 M3X8	-	Four
8	SET SCREW M3X5	-	One
9	SCREW, SST FLAT HEAD M3X10	-	Three




Fig. 10-2

10.3.3 Corona Needle ASSY

No.	Part Name	Part No.	Remarks
1	NEEDLE UNIT ASSY	S225-14290-41	Including needle and support arm
2	NEEDLE ASSY D	S225-15877-92	Needle only
3	HIGH VOLTAGE CABLE	S225-15888-41	Used in common with APCI
4	APCI/DUIS SOCKET ASSY	S225-14232-41	Including high-voltage cable




Fig. 10-3

10.3.4 IF Flange ASSY

No.	Part Name	Part No.	Remarks
1	HEATER FLANGE ASSY	S225-15486-41	
2	IF FLANGE	S225-15477-41	




Fig. 10-4

10.3.4.1 Heater Flange ASSY

No.	Part Name	Part No.	Remarks
1	SAMPLING CONE	S225-15487	
2	BLOCK, HEATER	S225-15488	
3	INSULATOR BLOCK	S225-15489	
4	HEATER ASSY	S225-15716-41	
5	INSULATOR, FLANGE	S225-15491	
6	BUSH, HEATER	S225-15498	
7	O-RING, 4D S100	S036-19004-53	
8	O-RING, 4D-S18	S036-19004-14	
9	O-RING, 4D S14	S036-19004-11	
10	O-RING, 4D-S10	S036-19004-07	
11	SCREW, SST FLAT HEAD M3X5	S225-14287-41	Two
12	BOLT M6X35	S225-15493-91	Two
13	BUSH FASTENER	S225-15494	Two
14	WASHER	S023-65106-01	Two




Fig. 10-5

10.3.4.2 IF Flange

No.	Part Name	Part No.	Remarks
1	DL PIPE ASSY	S225-15718-91	With O-rings (2 locations)
2	ORIFICE PART	-	
3	O-RING, 4D P4	S036-11202	




Fig. 10-6

■ DL (Desolvation Line)

No.	Part Name	Part No.	Remarks
1	O-RING, 4D P10	S036-11208	
2	O-RING, 4D-S6	S036-19004-03	

10




Fig. 10-7

■ Orifice unit

No.	Part Name	Part No.	Remarks
1	ORIFICE	S225-15479	
2	O-RING, 4D S8	S036-19004-05	
3	O-RING, 4D S14	S036-19004-11	
4	INSULATOR FLANGE	S225-15481	
5	PLATE, ORIFICE	S225-15710	
6	SCREW, SST PAN HEAD M2.5X3	-	Two




Fig. 10-8

10.3.5 Spray Unit

No.	Part Name	Part No.	Remarks
1	O-RING, 4D G145	S036-12527	
2	DOOR ASSY	S225-12100-41	
3	MESH	-	
4	DL PLUG	-	DL cover






Fig. 10-9

10.3.6 Lens System

No.	Part Name	Part No.	Remarks
1	Qarray ASSY	S225-12030-91	Qarray + skimmer, complete
2	SPACER	S225-12215-91	
3	OP1 ASSY	S225-12040-91	Multipole 1 + entrance lens, complete (for the LCMS-8030)
	QP1 ASSY	S225-13646-41	Multipole 1 + entrance lens, complete (for the LCMS-8040)
4	OP2 ASSY	S225-12050-91	Multipole 2 + entrance lens, complete (for the LCMS-8030)
	QP2 ASSY	S225-13648-41	Multipole 2 + entrance lens, complete (for the LCMS-8040)
5	O-RING, AS568A-253 4D	S036-15552-53	
6	HOOK	S225-12029	
7	KNOB, KNURL A-1176-17	S037-02829-12	




Fig. 10-10

10.3.6.1 Qarray ASSY: S225-12030-91

No.	Part Name	Part No.	Remarks
1	SKIMMER	S225-12034-01	
2	SCREW, SST FLAT HEAD M3X12	S020-12107	Two




Fig. 10-11

10.3.6.2 OP1 ASSY/QP1 ASSY

■ OP1 ASSY: S225-12040-91 (for the LCMS-8030)

No.	Part Name	Part No.	Remarks
1	OP1 LENS	S225-12045-01	
2	SCREW, SST SEMS P4 M3X6	S020-46634	




Fig. 10-12

■ QP1 ASSY: S225-13646-41 (for the LCMS-8040)

No.	Part Name	Part No.	Remarks
1	QP1 LENS	S225-13630-01	
2	SCREW, SST SEMS P4 M3X6	S020-46634	




Fig. 10-13

10.3.6.3 OP2 ASSY/QP2 ASSY

■ OP2 ASSY: S225-12050-91 (for the LCMS-8030)

No.	Part Name	Part No.	Remarks
1	OP2/QP2 LENS	S225-12052-01	
2	SCREW, SST FLAT HEAD M3X16	S020-12109	




Fig. 10-14

■ QP2 ASSY: S225-13648-41 (for the LCMS-8040)

No.	Part Name	Part No.	Remarks
1	OP2/QP2 LENS	S225-12052-01	
2	SCREW, SST FLAT HEAD M3X16	S020-12109	




Fig. 10-15

10.3.7 Detection Unit

No.	Part Name	Part No.	Remarks
1	DETECTOR	S225-14168	Detector, complete
2	ELECTRON MULTIPLIER	S225-14168-01	
3	Cable 1	S225-12292-91	For high voltage
4	Cable 2	S225-12293-91	For signals
5	SCREW, SST PAN HEAD M3X6	-	One
6	BOLT, SST HEXSOCH SEMS P3 M3X8	-	Two
7	SCREW, SST SEMS P4 M3X10	-	One




Fig. 10-16

10.3.8 Vacuum System

No.	Part Name	Part No.	Remarks
1	RP OIL Ultragrade19 (4 L)	S017-30163-02	About 1.5 L used at each oil change
2	IG GAUGE	S225-09490-01	Ion gauge vacuum gauge
3	FILAMENT FOR PB1, ROHS	S225-20310-91	Pirani gauge vacuum gauge
4	TMP, split flow 310	S225-14179-01	Triple inlet turbo molecular pump body + power supply
5	ROTARY PUMP, E2M28	S225-09309-02	For 230 V
6	CABLE, RP ASSY	S225-17224-41	Power supply cable for the rotary pump
7	LEAK VALVE ASSY	S225-03540-94	
8	HOSE TG-32	S018-31555-06	Use a 1.5 m RP intake hose (purchase unit: m)
9	HOSE, CHEMIFLEX 15MM	S016-31697-02	Use a 0.5 m TMP exhaust hose (purchase unit: m)
10	SLEEVE, PTFE 15X22 CL	S018-31511	RP exhaust hose (purchase unit: m)
11	HOSE BAND, GEAR S 50	S037-61064	For spring hose #32
12	HOSE BAND, WIRE 24	S037-61023	For spring hose #15
13	HOSE BAND, WIRE SY-22	S037-61002	For vinyl hose 15 × 1.5
14	CENTER RING, KF10	S035-06004-21	
15	CENTER RING, KF25	S035-06004-24	
16	CLAMP, KF16	S035-06004-01	
17	CLAMP, KF25C	S035-06004-02	
18	FLANGE PIPE, TQ-VAC	S225-12211-91	
19	ELBOW, KF25-#15	S225-03535-91	
20	FLANGE	S202-55374	For RP intake
21	O-RING, 4D P14	S036-11215	Drain seal
22	O-RING, 4D G50	S036-12506	Seal for the ionization probe
23	O-RING, 4D G100	S036-12517	Seal for the probe holder
24	O-RING, 4D P26	S036-11228	Seal for DUIS
25	O-RING, 4D P15	S036-11216	Pirani, hermetic seal
26	O-RING, 4D P40	S036-11243	For IG, feed-through
27	O-RING, AS568A-253 4D	S036-15552-53	OP cover
28	O-RING, AS568A-272 4D	S036-15552-72	Q3 cover
29	O-RING, AS568A-278 4D	S036-15552-78	Q1 cover
30	O-RING, AS568A-341 4D	S036-15553-41	Front and back housing
31	O-RING, 4D P18	S036-11218	For feed-through
32	O-RING, 4D P5	S036-11203	Feed-through conversion




Fig. 10-17

10.3.8.1 Ion Gauge Vacuum Gauge

No.	Part Name	Part No.	Remarks
1	IG GAUGE	S225-09490-01	
2	SACK NUT	S261-00210	
4	GASKET PRESS	S261-00209	
3	PACKING 4TX14X22	S261-00207-02	2 pc
5	FLANGE	S225-10150-01	
6	SCREW, SST SEMS P3 M4X12	-	Three
7	O-RING, 4D P40	S036-11243	For IG
8	BAFFLE	S225-15534	




Fig. 10-18

10.3.8.2 Pirani Gauge Vacuum Gauge

No.	Part Name	Part No.	Remarks
1	FILAMENT FOR PB1, ROHS	S225-20310-91	
2	O-RING, 4D P15	S036-11216	For the Pirani gauge
3	SCREW, SST SEMS P3 M4X8	-	




Fig. 10-19

10.3.9 Waste Tube

No.	Part Name	Part No.	Remarks
1	SILICON RUBBER TUBE, 7X10NL	S016-31350-19	1 m drain tube from the leak tray
2	PVC TUBE, R3603 1/2X3/4X1/8	S016-31414	1 m drain tube from the spray unit, for use outside the instrument
3	TUBE, SE-200 1/2-3/4	S016-37619-02	320 mm drain tube from the spray unit, for use inside the instrument
4	SLEEVE, PTFE 12X3	S018-31510	5 m
5	RUBBER CAP	S225-06482-92	




Fig. 10-20

10.3.10 Standard Sample Introduction Unit

No.	Part Name	Part No.	Remarks
1	FRIT PEEK	S228-48607-91	One
2	CAPILLARY ASSY	S225-15848-91	
3	Standard sample bottle	S038-00512-01	
4	SPACER, FKM	S225-15697-01	
5	Bottle cap	S225-15868-91	
6	Standard Sample	S225-14122-01	Standard sample solution 200 mL
7	BRACKET, SI	S225-15560	
8	SCREW, SST FLAT HEAD M3X6	-	
9	HALF UNION, GWJS6-M5	S035-65415-06	
10	HOUSING, LINEFILTER	S228-46358	
11	MALE NUT, PEEK	S228-18565-84	5 included
12	PEEK TUBE 1.6X0.25	S228-32999-03	3 m (cut to 95 mm for use)
13	PEEK TUBE 1.6X0.25	S228-32999-03	3 m (cut to 400 mm for use)
14	TUBE, 0624304	S225-15873-91	Resistance tube
15	MAIL CONNECTER, PEEK	S228-25014	




Fig. 10-21

10

10.3.11 Gas Controller

No.	Part Name	Part No.	Remarks
1	Filter, 2300B-SS-1/8-2U	-	
2	VALVE ASSY (NEB)	S040-51226-01	
3	VALVE ASSY (DRY)	S225-14292-41	
4	SUS316L 1.6X1.0X40	S228-49000-40	
5	Valve, SI	S225-15688-91	
6	VALVE, RELIEF RAB2V-150	S040-27013-01	
7	REGULATOR, AR10-M5BG	S040-72549-51	
8	VALVE ASSY SAGINOMIYA	S221-48813-91	
9	FILTER,PTFE 20 MICRON	S221-18154	
10	QUARTZ TUBE, TSP040375	S016-70402-01	
11	VALVE, RELIEF	S225-15688-42	
12	LATERAL, UNION KQ2U06-00	S035-60693-29	
13	FEP TUBE 6 mm, GAS	S225-15846-93	
14	FEP TUBE 1/16, NEB	S225-14255-41	
15	FEP TUBE 4 mm, DRY	S225-14255-42	
16	FEP TUBE 1/16, DRY	S225-14255-43	
17	POLYURETHANE TUBE, U2-4-6X4BK	S016-46021	
18	TEE, RUN KQ2Y04-M5-ROHS	S035-60690-05	
19	TEE, RUN KQ2Y06-M5-ROHS	S035-60690-08	
20	TUBE, M5-MF 80MM	S221-41410-93	
21	3-WAY JOINT ASSY	S221-25211-91	
22	TUBE, M5-MF 260MM	S221-41410-94	
23	TUBE, M5-MF 200MM	S221-41410-92	




Fig. 10-22

10.3.12 Accessory Parts

No.	Part Name	Part No.	Remarks
1	SPANNER, HEX 1.5 ROHS	-	
2	SPANNER, HEX 3 ROHS	-	
3	SPANNER, HEX 5 ROHS	-	
4	SPANNER, HEX 6 ROHS	-	
5	SPANNER, HX7-8-R	-	
6	SPANNER, HX10-12-R	-	
7	SPANNER, DOUBLE OPEN END17X19-R	-	
8	SCREWDRIVER, TORQUE #2 100MM	-	
9	CABLE, LCMS EVENT	S225-17126-41	
10	CONNECTER, 6 P	-	
11	CONNECTER, 8 P	-	
12	PEEK CAPILLARY 5999B	-	
13	DL PLUG	-	
14	JIG PULL	-	
15	MALE NUT, PEEK	S228-18565-84	
16	SPACER, FKM	S225-15697-01	
17	MAGNIFIER, 1962	-	
18	CLAMP, DKN-5GSP	-	
19	GUARD	-	
20	ABRASIVE CLOTH	-	
21	COUPLING, 1.6C 316L	-	




Fig. 10-23

10.3.13 Startup Kit (Option)

No.	Part Name	Part No.	Remarks
1	FERRULE 1.6F	S228-33513-91	PEEK ferrule (total of 3)
2	DL PIPE ASSY	S225-15718-91	
3	PUMP OIL, H11025013	S017-30163-02	Rotary pump oil, 4 L
4	MALE NUT, PEEK	S228-18565-84	Male nut, 1.6 MN (total of 5)
5	FERRULE, 1.6F-T	S228-16007-84	PTFE ferrule (total of 5)
6	PEEK CAPILLARY 5999B	-	ID0.13 mm PEEK tube, 3 m
7	ETFE TUBE 1.6X0.5	S228-18495-04	
8	FITTING	S228-32651-41	Nut, ferrule (5 each)
9	CAPILLARY ASSY	S225-14948-91	
10	SI BOTTLE LINE FILTER	S228-48607-91	
11	0.1 SUS PIPE	S228-49120-00	ID0.1 mm SUS pipe, 2 m
12	PEEK TUBE CUTTER	S228-32930-01	
13	FILTER	S042-60935-14	
14	GROVE, LATEX	S086-72599-01	
15	RESISTANCE TUBE, GLP CLEAN	S228-32722-94	




Fig. 10-24

10.3.14 RP Oil Return Kit (Option): S225-05990-92

No.	Part Name	Part No.	Remarks
1	OIL RETURN KIT, E2M28	-	
2	OIL MIST FILTER, EMF20	S042-00124-33	
3	FLANGE, KF25-#15 STRATE	-	
4	CENTER RING, KF25ANCR	S035-06004-13	
5	CLAMP, KF25C	S035-06004-02	
6	SPANNER, HEX 8 ROHS	-	




Fig. 10-25

Index

A

Ammonium acetate	213
Ammonium ions	215
Analysis unit	28, 211
Anoxia symptoms	193
APCI heater	94
APCI needle	136
APCI pipe	129
APCI probe	20, 128
APCI set	8
APCI/DUIS socket	9
Ar gas supply port	16
Argon gas	193
Atmospheric pressure ionization	18
Auto tuning	80
Automatic restart mode	116

C

Cable tie	6
Capillary	122
Carrier gas tubing	196
CBM-20A	192
CID gas	16, 195
Collision cell	28, 30, 36, 99
Collision gas control valve	36
Collision gas stop valve	36
Collision-induced dissociation	101
Conversion dynode	31, 32
Corona discharge	20
Corona needle	136

D

Deprotonated molecules	215
Detector	17, 31, 212
Detector adjustment	84
Dimer ions	215
Distribution panel	105
DL	23, 94, 163
Drying gas	16, 22, 96
Dual ion source (DUIS)	72
DUIS (dual ion source)	21

DUIS analysis	72
DUIS high-voltage cable connector	74
DUIS needle	140
DUIS set	9

E

Electron multiplier	31
Entrance lens	25, 27
ESI analysis	55
ESI capillary tube	55
ESI coupling	124
ESI high-voltage cable connector	21
ESI probe	19, 55, 58, 75, 122

F

Fan	181
FCV-20AH	16
FCV-20AH mounting kit	12
Ferrule	134
Filter	157
Flow rate monitor	96

G

GAS	38
Gas	96, 193
Gas ballast valve	175
Gas joint	138
Gas pipe adapter kit	196

H

Half union	138
Heated block	22, 94
HEATER	38
Heater	94
Heater cable connector	20
Heater flange	164
Heater unit	132
High voltage	93
High-pressure gas cylinder	194
High-pressure gradient system	216

H

High-voltage cable connector 74

I

I/O signal connector 16

IF flange 165

Inch size 196

Indication lamps 38

Ion gauge 33

Ion guide mode 99

K

Knee region 32

L

LC 184, 192

Leak sensor 37, 180

Leak tray 180

Leak valve 33, 114, 115

LEDs 38

Lens 25, 142

Lens voltage 210

Line filter 34, 36

Loupe 4

M

Manual tuning 201

Mass calibration 85

Mass chromatogram 100

Mass spectrum 99

MC 100

Mist exhaust hose 189

MRM 103

MS detector 31

MS scan mode 99

MS/MS analysis mode 101

Multipole 25, 27

N

N₂ gas supply port 16

Nebulizer gas 16, 21, 96

Nebulizer gas tube 178

Needle alignment tool 73

Needle unit 136, 140

Negative ionization mode 215

Neutral fragments 102

Neutral loss scan mode 102

Nitrogen gas 193

O

Oil filler plug 172

Oil mist filter 190

Oil return tube 190

One-touch joint for gas 57, 195

Orifice 22, 24, 166

P

Peak profile measurement 85

PEG 215

Pirani gauge 33, 90

Plateau region 32

Plug DL 4

Polyethylene glycol 213

Polypropylene glycol 213

Position adjustment knob 56

Positive ionization mode 215

Post-rod 28

Precursor ion 101

Pre-rod 28

Probe holder 145, 149

Product ion 102

Prominence system 216

Pseudo-potential 30

Q

Q1 post-rod 2, 17

Q1 pre-rod 17

Q1 quadrupole rod 2, 17

Q2 collision cell 2, 17

Q3 pre-rod 2, 17

Q3 quadrupole rod 2, 17

Qarray 2, 17, 25, 26

Quadrupole rod 28, 93

R

Raffinose	213, 215
Release bush	57
Resistance tube	34, 36, 156
Resolution adjustment.....	84
Resolution mode	99
Rotary pump.....	171
RP exhaust hose set.....	6
RP oil return kit.....	10

S

Sampling cone	22
Scan speed	100
Secondary electron	31
Secondary electron multiplier.....	31
Secondary electrons	32
Selective reaction monitoring	103
Sensitivity adjustment	84
Shimadzu M type joint.....	195
Shimadzu M-shaped joint.....	16
SIM mode.....	100
Skimmer	2, 17, 25, 26
Spray unit.....	159
Standard sample bottle	14, 34
Standard sample introduction unit.....	34
Start signals	16
Startup kit.....	7, 11
Stopping the instrument	48
System controller	192

T

TIC chromatogram (total ion current chromatogram).....	99, 100
Triple inlet turbo molecular pump.....	33
Tuning file.....	207

U

USB.....	15
User Management.....	220

V

Vacuum system.....	33
Vinyl tube	6

W

Waste drain	17, 37
Waste tube	179

This page is intentionally left blank.